

**Materials List for:**

# Measuring the 50% Haemolytic Complement (CH<sub>50</sub>) Activity of Serum

Maurizio Costabile<sup>1</sup>

<sup>1</sup>School of Pharmacy and Medical Sciences, University of South Australia

Correspondence to: Maurizio Costabile at [maurizio.costabile@unisa.edu.au](mailto:maurizio.costabile@unisa.edu.au)

URL: <https://www.jove.com/video/1923>

DOI: [doi:10.3791/1923](https://doi.org/10.3791/1923)

## Materials

Name	Company	Catalog Number	Comments
Sheep red blood cells	CSL	2490201	Use at 1% final
Serum			Human serum
Rabbit Anti-Sheep Haemolytic serum (RBCs), Unconjugated	AbD Serotec	C12HSB	
96 well flat bottom plate	Sarstedt Ltd	83.1839	
Veronal buffered saline			Use at 1x final
Waterbath			
Plate reader			
37°C room/incubator			