

Video Article

Flash Freezing and Cryosectioning E12.5 Mouse Brain

D. Spencer Currle¹, Edwin S. Monuki¹

¹Department of Developmental and Cell Biology, University of California, Irvine (UCI)

Correspondence to: D. Spencer Currle at Spencer.Currle@JoVE.com

URL: https://www.jove.com/video/198

DOI: doi:10.3791/198

Keywords: Neuroscience, issue 4, mouse, brain, sectioning

Date Published: 5/28/2007

Citation: Currle, D.S., Monuki, E.S. Flash Freezing and Cryosectioning E12.5 Mouse Brain. J. Vis. Exp. (4), e198, doi:10.3791/198 (2007).

Abstract

Video Link

The video component of this article can be found at https://www.jove.com/video/198/

Protocol

- 1. Fix tissue in 4% paraformaldehyde in PBS for desired time.
- Sucrose infuse tissue (cryoprotection)
 - 1. Make 30% sucrose solution in PBS w/v in 2059 tube.
 - 2. Rinse tissue 3x in PBS (~ 5 min with rocking).
 - 3. Place tissue in 30% sucrose solution. Tissue will not sink.
 - 4. Place the tissue in 4°C overnight, or until it has sunk.
- 3. Label appropriate size cryomold with information and orientation.
- 4. Fill cryomold with O.C.T. (avoid bubbles).
- 5. Transfer tissue to O.C.T. bath and coat it with O.C.T.6. Transfer tissue to O.C.T. in cryomold.
- 7. Orient the tissue under microscope.
- 8. Pour liquid nitrogen into plastic Petri dish.
- 9. Quickly and carefully lower the tissue in cryomold into the nitrogen. (do not submerge the top of the cryomold.)
- 10. When the O.C.T. is solid white, place the frozen tissue into -80°C freezer for storage.
- 11. Equilibrate tissue to ~20°C for at least 30 min. prior to sectioning.