

Materials List for:

Streamlined Purification of Plasmid DNA From Prokaryotic Cultures

Laura Pueschel¹, Hongshan Li¹, Matthew Hymes¹

¹Pall Life Sciences

Correspondence to: Matthew Hymes at matthew_hymes@pall.com

URL: <https://www.jove.com/video/2407>

DOI: [doi:10.3791/2407](https://doi.org/10.3791/2407)

Materials

Name	Company	Catalog Number	Comments
AcroPrep Advance Filter plate	Pall Corporation	PN 8175 (1 mL) or PN 8075 (350 μ l)	3 μ m glass fiber 0.2 μ m Supor® membrane for clarification
AcroPrep Advance Filter Plate	Pall Corporation	PN 8132	for DNA Binding
Vacuum manifold	Pall Corporation	PN 5017	
Vacuum pressure pump	Pall Corporation	PN 13157	