

Materials List for:

Dissection of a Mouse Eye for a Whole Mount of the Retinal Pigment Epithelium

Alison Claybon¹, Alexander J. R. Bishop¹

¹UT Health Science Center at San Antonio, Greehey Children's Cancer Research Institute and Department of Cellular and Structural Biology

Correspondence to: Alexander J. R. Bishop at bishopa@uthscsa.edu

URL: <https://www.jove.com/video/2563>

DOI: [doi:10.3791/2563](https://doi.org/10.3791/2563)

Materials

Name	Company	Catalog Number	Comments
straight forceps	Roboz Surgical Instruments Co.	RS-4903	tip: .08 x .04 mm material: INOX
45° forceps	Roboz Surgical Instruments Co.	RS-5005	tip: .05 x.01 mm material: INOX
15° "up" forceps	Roboz Surgical Instruments Co.	RS-5045A	tip: .1 x.06 mm material: INOX
spring scissors	Roboz Surgical Instruments Co.	RS-5604	comb. tip width 0.3mm cutting edge length 3mm material: stainless steel
binocular dissecting microscope	Carl Zeiss, Inc.	Discovery V.8	use reflected light source
phalloidin	Invitrogen	A22283	Alexa Fluor 546