

Materials List for:

Extracting DNA from the Gut microbes of the Termite (Zootermopsis Angusticollis)

Eric Matson¹, Jared Leadbetter¹

¹Department of Environmental Science and Engineering, California Institute of Technology - Caltech

Correspondence to: Jared Leadbetter at jleadbetter@caltech.edu

URL: <http://www.jove.com/video/389>

DOI: [doi:10.3791/389](https://doi.org/10.3791/389)

Materials

Name	Type	Company	Catalog Number	Comments
PVPP/SDS/phenol buffer	Buffer			homogenization buffer
DNeasy	kit	Qiagen		
TE	Buffer			1x buffer: 1 mM Tris-HCl, 0.1 mM EDTA, pH 8.0
zirconia/silica beads				0.1 mm
PVPP				1x TE buffer containing 1% w/v polyvinylpyrrolidone (PVPP)
Zootermopsis Angusticollis	Animal			Termites