

Video Article

Visualization of Cortex Organization and Dynamics in Microorganisms, using Total Internal Reflection Fluorescence Microscopy

Felix Spira¹, Julia Dominguez-Escobar¹, Nikola Müller^{1,2}, Roland Wedlich-Söldner¹

¹AG Cellular Dynamics and Cell Patterning, Max Planck Institute of Biochemistry

²Helmholtz Zentrum München

Correspondence to: Roland Wedlich-Söldner at wedlich@biochem.mpg.de

URL: <http://www.jove.com/video/3982>

DOI: [doi:10.3791/3982](https://doi.org/10.3791/3982)

Keywords: Bioengineering, Issue 63, TIRF, imaging, yeast, bacteria, microscopy, cell cortex

Date Published: 5/1/2012

Citation: Spira, F., Dominguez-Escobar, J., Müller, N., Wedlich-Söldner, R. Visualization of Cortex Organization and Dynamics in Microorganisms, using Total Internal Reflection Fluorescence Microscopy. *J. Vis. Exp.* (63), e3982, doi:10.3791/3982 (2012).

Abstract

TIRF microscopy has emerged as a powerful imaging technology to study spatio-temporal dynamics of fluorescent molecules in vitro and in living cells¹. The optical phenomenon of total internal reflection occurs when light passes from a medium with high refractive index into a medium with low refractive index at an angle larger than a characteristic critical angle (i.e. closer to being parallel with the boundary). Although all light is reflected back under such conditions, an evanescent wave is created that propagates across and along the boundary, which decays exponentially with distance, and only penetrates sample areas that are 100-200 nm near the interface². In addition to providing superior axial resolution, the reduced excitation of out of focus fluorophores creates a very high signal to noise ratios and minimizes damaging effects of photobleaching^{2,3}. Being a widefield technique, TIRF also allows faster image acquisition than most scanning based confocal setups.

At first glance, the low penetration depth of TIRF seems to be incompatible with imaging of bacterial and fungal cells, which are often surrounded by thick cell walls. On the contrary, we have found that the cell walls of yeast and bacterial cells actually improve the usability of TIRF and increase the range of observable structures⁴⁻⁸. Many cellular processes can therefore be directly accessed by TIRF in small, single-cell microorganisms, which often offer powerful genetic manipulation techniques. This allows us to perform *in vivo* biochemistry experiments, where kinetics of protein interactions and activities can be directly assessed in living cells.

We describe here the individual steps required to obtain high quality TIRF images for *Saccharomyces cerevisiae* or *Bacillus subtilis* cells. We point out various problems that can affect TIRF visualization of fluorescent probes in cells and illustrate the procedure with several application examples. Finally, we demonstrate how TIRF images can be further improved using established image restoration techniques.

Video Link

The video component of this article can be found at <http://www.jove.com/video/3982/>

Protocol

1. Sample Preparation

a. Preparing cover slips

Coverslips should be cleaned from dust particles as TIRF is sensitive to background signals arising from the coverslip (**Fig. 1A, Movie 1**).

1. Using forceps place cover slips in a ceramic holder with lid.
2. Fill glass container with 1 M NaOH (can be reused multiple times).
3. Incubate for 2 h under slow continuous rotation (orbital shaker). Excessive incubation (> 8 h) will create opaque cover slips.
4. Wash twice with distilled water for 5 min under slow continuous rotation.
5. Store in ceramic holder covered in 100% ethanol.

b. Preparing *Bacillus subtilis* cells

1. Dilute to an OD₆₀₀ of 0.01 - 0.05 in 5 ml of appropriate growth media and grow to exponential phase (OD₆₀₀ of 0.3-0.7).
2. Prepare 1.25% agarose in water or medium. Use synthetic media to reduce background fluorescence. Dissolve powder in a 1.5 ml plastic tube at 95°C and store in heating block at 72°C.
3. Add a small drop of agarose to the middle of a slide and press into a flat pad with a second slide. Carefully separate slides after at least 2 min or before use.
4. Spin down 500 µl cells in micro centrifuge at 12.000 rpm for 1 min.
5. Discard supernatant, and resuspend pellet in 50 µl medium.

6. Add 2-4 μ l cells to the centre of the agarose pad.
7. Take a cleaned cover slip from the ethanol filled container with forceps, dry with pressurized air and carefully place on sample.
8. Let cells settle for at least 2 min prior to microscopy.
9. Seal edges of coverslip with VALAP (Vaseline, Lanolin and Parafin, mixed at equal weight under heat, apply with small brush or spatula) for long term imaging.

c. Preparing *Saccharomyces cerevisiae* cells

1. Dilute pre-culture and grow in appropriate media for at least 5 h to exponential phase (OD_{600} 0.2-0.8).
2. Take a coverslip from the container with forceps and dry with pressurized air.
3. Spread 5 μ l Concanavalin A solution (ConA, 2 μ g/ml in distilled water) on cover slip with a pipette and dry with pressurized air. ConA binds to carbohydrates of the yeast cell wall and immobilizes yeast cells.
4. Transfer 4.5 μ l yeast cell suspension to one side of the ConA-coated coverslip. Carefully place cover slip on microscope slide. Start with one edge and let drop slowly to avoid formation of air bubbles.
5. Let cells settle for at least 2 min prior to microscopy.
6. Seal edges of coverslip with VALAP (see above) for long term imaging.

2. Image Acquisition

a. Microscope setup

We performed all experiments on a customized TIRF setup based on a fully automated iMIC stand (Till Photonics) with an Olympus 100x 1.45 NA objective. 75 mW DPSS lasers at 488 nm (Coherent Sapphire) and 561 nm (Cobolt Jive) were used as light sources. Lasers were selected through an AOTF and directed through a broadband fiber to the iMIC. A galvanometer-driven two-axis scanner head was used to adjust incidence angles or Fluorescence Recovery after Photobleaching (FRAP) position. An additional galvanometer was used to switch between epifluorescence, FRAP and TIRF. Images were collected with an Andor iXON DU-897 EMCCD camera and a 2x magnification lens in front of the camera. Acquisition was controlled by the Live Acquisition (Till Photonics) software package.

The following points should be considered:

1. Programmable control over illumination source (laser line) and TIRF angle is essential for reproducible results and especially for multi-color TIRF, where incidence angles need to be adjusted for each excitation wave length to ensure equal penetration of evanescent waves. In our setup TIRF angle can be adjusted instantaneously via two galvanometers, which are directly controlled from the Live Acquisition software package.
2. Our customized setup additionally allows rapid switching between epifluorescence, TIRF and FRAP via a third galvanometer, allowing simple measurement of localization kinetics for objects traced in TIRF.
3. Objective-type TIRF requires high numerical apertures of $NA > 1.4$. We used an Olympus 100x NA 1.45 objective.
4. EMCCD cameras work very well with TIRF due to their high signal to noise ratio. To allow imaging with optimal sensitivity we used a 512x512 pixel back illuminated EMCCD camera with a pixel size of 16 μ m. To maintain maximum resolution we placed a 2x magnification lens in front of the camera.
5. As TIRF is a widefield technique, intensity of laser beams is spread out. Strong lasers therefore often improve imaging. We used 75 mW diode pumped solid state (DPSS) lasers at 488 nm and 561 nm.
6. To keep cells under optimal growth conditions during imaging, use appropriate environmental controls. We used a custom made heating chamber (Fig. 1B) to study yeast temperature sensitive mutants. For medium exchange during acquisition we either used simple flow chambers or glass-bottom culture dishes that can be easily accessed from above on the inverted microscope.

b. Adjusting parameters

1. Identify position of cells using bright field illumination. Red LEDs are particularly good as they do not induce significant photo damage.
2. Switch to laser illumination and select appropriate lasers and filter combinations to excite the fluorophores of choice. Make sure that the TIRF angle is subcritical, or you will not be able to detect any signal arising from GFP fusion proteins. If in doubt set TIRF angle to 0° to obtain a perpendicular laser beam. Find the top section of the cell (side of the cell facing the coverslip). Alter the incident angle of the laser beam until signals disappear then slowly increase angle until the point where fluorescence at the cell surface reappears (**Movie 2**). Adjust z-focus again for optimal position at the surface. Acceptable TIRF angles create no blurred halo at the edge of the cell (**Fig. 1C**). If significant photo bleaching occurs during setup, save the TIRF angle and move to an unbleached stage position before starting acquisition.
3. Adjust laser intensities and camera gain to maximize dynamic range (high signal to noise ratio without saturating pixels). Typically EMCCD cameras are optimal to detect very low signals at high signal to noise ratios. For stronger signals regular CCD cameras produce less noise and offer a wider dynamic range.
4. TIRF microscopy is very sensitive to small height differences of the coverslip, resulting in only few cells, which can be imaged at the same time (**Fig. 1D**). For small cells the acquisition region (region of interest, ROI) can therefore often be reduced to a subregion containing the object of choice. This also reduces the time needed to transfer data from the chip to the hard drive, often the rate-limiting step at high frame rates.
5. For double color TIRF experiments make sure to minimize bleed through of fluorophores. For RFP/GFP pairs use separate emission filters as RFP is weakly excited by 488 nm light (**Fig. 2A**). Most filters are not perfectly aligned with each other. Filter offset can be corrected by using fluorescent beads mixed in with cells. To achieve comparable penetration depth, adjust TIRF angles separately for each laser line. A typical workflow is illustrated in **Fig. 2B**.

6. A critical parameter influencing image quality is laser alignment and a clean objective. To align the laser, first focus on the z-position used for TIRF imaging. Then remove the sample and clean the objective of all oil (residual oil will lead to diffraction of the laser light and speckles in the beam profile). Project the laser on the ceiling in TIRF mode and calibrate the 0° position such that the laser is in a straight line with the objective (use reference spot, always wear laser safety goggles and avoid all reflections in the beam path). Focus the laser spot to a minimal size. Our setup provides a moveable lens positioned between galvos and microscope to facilitate rapid calibration. After optimal adjustment, the laser profile should form a well define spot with roughly round shape. Perform calibration before every session to obtain optimal image quality.

3. Representative Results

- a. **Distribution of actin homologues and cell wall enzymes in *Bacillus subtilis*.** We imaged cells expressing GFP fusions of the actin homologue Mbl or the transpeptidase PbpH by TIRF and regular epifluorescence. The penetration depth in TIRF is clearly reduced and limited to the surface facing the coverslip (**Fig. 3A, B**). The weakly expressed transpeptidase PbpH localizes to cortical patches that are hardly visible by epifluorescence but can be clearly distinguished by TIRF (**Fig. 3B**). The reduced penetration can be best observed for PbpHGFP signal at septa, which appear as lines in epifluorescence but as dots in TIRF (**Fig. 3B, arrow**).
- b. **Protein domains in the yeast plasma membrane.** Examples for network-like patterns formed by plasma membrane proteins in yeast are already given in **Fig. 1C and 2**. In addition to the improved contrast in TIRF that allows clear distinction of these network-like patterns, weak signals can sometimes be detected exclusively by TIRF microscopy. As example, we imaged a GFP fusion of the TORC2 complex component Bit61 from *Saccharomyces cerevisiae*. This protein is weakly expressed and forms small cortical patches with only few protein copies per patch⁹. In epifluorescence only few patches are visible above the strong background, while patches can be imaged with very good signal to noise ratio in TIRF (**Fig. 3C**). Therefore TIRF is particularly well suited to study weakly fluorescent structures and proteins at the cell surface.

The high axial resolution provided by TIRF already reduces background fluorescence. Additional improvement of image contrast can be obtained through various image processing steps. A simple standard technique is local background equalization where a blurred image is subtracted from the original. Blurring can be performed with a variety of low pass filters including median or Gaussian filters (**Fig. 3D**). Local background subtraction removes noise and sharpens high intensity structures. However, this often comes at the price of a loss of fine structures and exaggerated amplification of high intensity regions (**Fig. 3D, arrow**). A superior procedure is 2D deconvolution, which can be performed with free or commercially available software packages such as ImageJ, Matlab or Huygens (Scientific Volume Imaging b.v.). TIRF images provide very good axial resolution and therefore a nearly two-dimensional point spread function (PSF). We experimentally determined this PSF by acquiring images of 40 nm fluorescent beads with the same settings that were used for the actual imaging (objective, camera, excitation and emission wave length). The experimental PSF was then used for deconvolution with a classical maximum likelihood algorithm in Huygens. Time lapse movies of GFPRas2 or of Fet3GFP and Pma1RFP (**Movies 3 and 4**) demonstrate the increase of contrast after deconvolution. In double color images it is critical to maximize contrast, to be able to quantify colocalization values (**Movie 4**).

- c. **Analysis of protein turnover at the cell cortex.** TIRF and FRAP offer a powerful combination for the study of dynamics of cortical proteins. However, to allow the efficient use of these techniques, rapid control of laser angles and positions is required. In TIRF, a wide field technique, the laser needs to be focused on the back focal plane and needs to follow a shallow incidence angle to achieve the desired reflection. In contrast, FRAP requires the laser to be focused directly on the sample and positioned with high spatial accuracy to allow bleaching of defined structures or regions. The setup from Till photonics we used in our experiments provides very fast and intuitive control over these parameters. Two galvanometer-controlled mirrors are used to adjust the FRAP position in x,y or to set the incidence angle in TIRF. A third mirror is used to switch between separate beam paths for TIRF and FRAP illumination (the laser is widened by additional lenses in the TIRF path). All settings are controlled by the LA (live acquisition) software and calibrations can be rapidly performed for each experiment to achieve optimal conditions. This software also allows definition of bleach positions during image acquisition, which is essential for FRAP experiments on motile structures.

Using this feature we could exclude treadmilling of Mbl filaments as source for the observed motility of Mbl-containing patches (**Fig. 3E**). This observation motivated our search for alternative motility mechanisms, and finally resulted in our identification of a completely unexpected type of motility, where intracellular protein complexes (including Mbl) are moved through the activity of cell wall synthesizing enzymes residing on the outside of the cell⁴. In a similar fashion, we were able to observe the slow rearrangements of yeast plasma membrane proteins, which distribute into network-like domains covering the whole cell surface (**Fig. 3F**).

These examples illustrate how the dynamics and patterning of cortical proteins can be directly assessed by combining TIRF and FRAP techniques.

4. Representative Results




Figure 1. Adjusting parameters. **A.** Dirty vs. cleaned coverslip. Background signal arising from dust on the coverslip interfere with a potential GFP/RFP signal, cleaned coverslip has less background. **B.** Custom-made electrical heating controller showing set temperature (a), temperature at probe (b) and sample holder (c). **C.** Yeast protein Pma1GFP imaged with decreasing incidence angles (from top left to bottom right). Images show a sudden loss of signal at the critical angle and progressive decrease in structural information and signal to noise ratio. **D.** Uneven field of view. A dense suspension of fluorescent beads is imaged, showing that only part of the field of view is in focus. Scale bars in A, B: 2 μ m, in C: 10 μ m.




Figure 2. Double color imaging. **A.** Bleed through of Pil1RFP. Pil1RFP is excited with 488 nm and a 561 nm laser and fluorescence is recorded with a double band pass filter. Signal of Pil1RFP is weakly visible in the GFP channel. Scale bar 2 μ m **B.** Typical steps to avoid bleed through in double color imaging. After imaging the cells with separate filter sets, they are deconvolved and aligned (using fluorescent beads) before being merged for analysis.




Figure 3. Representative Results. **A.** Comparison of GFP-Mbl distribution in *B. subtilis* using TIRF or regular epifluorescence (Epi) microscopy. Blue: Cell boundary from bright field image. Images represent average projections of a time series to indicate the area covered by motile MreB patches monitored with different imaging modes. **B.** Improved axial resolution and contrast of TIRF images taken of the transpeptidase PbpHGFP in *B. subtilis*. Arrows indicate septal staining that appears as ring in epifluorescence and as patches in TIRF. Exposure time: 100 ms. **C.** Comparison between epifluorescence and TIRF illumination of the TOR complex component Bit61 in *S. cerevisiae*. Bit61 is very weakly expressed⁹. Exposure time: 250 ms. **D.** Comparison of different image processing methods. Series showing raw TIRF image of Pma1GFP, subtraction of Gaussian or Median blurred images, as well as deconvolution using the maximum likelihood algorithm in Huygens. **E.** TIRF-FRAP of a single GFP-Mbl patch (above asterisk) moving across a *B. subtilis* cell. The kymograph of the non-recovering patch and the intensity profile along the kymograph (dotted line) rule out treadmilling as source of motion. **F.** TIRF-FRAP of a yeast cell expressing Pma1GFP. Note the diffusive closing of the gap in the cortical staining pattern. Scale bars: 2 μ m. Time in seconds.

Movie 1. Movie showing *B. subtilis* cells expressing an RFP fusion protein on a dirty coverslip. Note that RFP signal is hardly distinguishable from background noise. Cycle time 100 ms. [Click here to view movie.](#)

Movie 2. Movie of Pma1GFP with different incidence angles. Angle is altered from subcritical to critical angle, until signal is lost. Note, at subcritical angle internal signal is visible, causing noisy images until at critical angle just the protein at the PM is visible. Cycle time 200 ms. Scale bar: 2 μ m. [Click here to view movie.](#)

Movie 3. Movie of yeast GFPRas2, showing alternating RAW and deconvolved TIRF images. GFPRas2 is a fast moving protein (T1/2 of 2s unpublished data), fast acquisition times are important to resolve the dynamic behavior. Cycle time 150 ms. Scale bar: 2 μ m. [Click here to view movie.](#)

Movie 4. Double color movie of yeast proteins Fet3GFP and Pma1RFP. First 10 frames represent RAW images and last frames are deconvolved. This movie shows the importance of image processing for colocalization analysis. The blurry RAW images become sharpened, making an analysis possible. Cycle time 5 s. Scale bar: 2 μ m. [Click here to view movie.](#)

Discussion

TIRF microscopy is the technique of choice to image peripheral proteins. The low penetration depth of the evanescent field minimizes of out of focus light, which leads to a very good signal to noise ratio and allows data acquisition with high frame rates, or imaging of very weakly

expressed proteins. The combination of high contrast and high frame rates makes TIRF microscopy a perfect tool for imaging spatio-temporal dynamics of cortically localized proteins. Interestingly the thick cell wall surrounding many microorganisms does not hamper imaging of peripheral proteins by TIRF. In fact, the actual generation of the evanescent very likely occurs at the interface between cell wall and plasma membrane^{5,10}. Reflection of light that exits the cell might even lead to propagation of light along the cell wall, which might explain the large surface area that could be imaged in yeast cells⁸.

For optimal quality of TIRF images it is crucial to have a well calibrated microscopy system. The laser should be focused and aligned before each imaging session. Coverslips should be cleaned (**Fig. 1A**) and the cells have to be properly immobilized. For double color imaging spectral bleed through has to be taken into account or eliminated by using filter sets filters separate emission range. This will obviously come at the price of slower acquisition times due to mechanical filter change. Fast filter wheels or galvanometer driven illumination sources can circumvent this delay if necessary. Alternatively, the critical fluorophore (RFP) can be used on the weaker expressed protein in a pair, minimizing the level of signal interference.

Despite of the high contrast, images taken with a TIRF microscope contain significant noise. To remove this noise and further increase image contrast images can be processed with various filtering steps. The method with the best results in our experience is 2D deconvolution. This requires the experimental measurement of the PSF for each sample and microscope condition. However, this can be performed fairly easily and quickly using fluorescent beads mixed into the respective sample. For two-color experiments these beads can additionally be used for channel alignment.

We have demonstrated the advantages of using TIRF microscopy for imaging of cortical proteins in microorganisms. Combination of TIRF and image processing allows acquisition of images with high frame rates (<50 ms) and contrast and also enables visualization of weakly expressed proteins such as Bit61. An additional advantage of performing TIRF imaging on microorganism is that turgor pressure in these cells leads to flat plasma membranes, and minimizes topological effects to TIRF signals.

Disclosures

No conflicts of interest declared.

Acknowledgements

This work was funded by the Max Planck society.

References

1. Axelrod, D., Thompson, N.L., & Burghardt, T.P. Total internal reflection fluorescent microscopy. *J. Microsc.* **129**, 19-28 (1983).
2. Axelrod, D. & Omann, G.M. Combinatorial microscopy. *Nat. Rev. Mol. Cell Biol.* **7**, 944-952, doi:10.1038/nrm2062 (2006).
3. Axelrod, D. Selective imaging of surface fluorescence with very high aperture microscope objectives. *J. Biomed. Opt.* **6**, 6-13, doi:10.1117/1.1335689 (2001).
4. Dominguez-Escobar, J., *et al.* Processive movement of MreB-associated cell wall biosynthetic complexes in bacteria. *Science*. **333**, 225-228, doi:10.1126/science.1203466 (2011).
5. Sparkes, I.A., *et al.* Bleach it, switch it, bounce it, pull it: using lasers to reveal plant cell dynamics. *J. Exp. Bot.* **62**, 1-7, doi:10.1093/jxb/erq351 (2011).
6. Uchida, M., *et al.* Total synthesis and absolute configuration of avenolide, extracellular factor in *Streptomyces avermitilis*. *J. Antibiot. (Tokyo)*, doi:10.1038/ja.2011.90 (2011).
7. Yu, H. & Wedlich-Soldner, R. Cortical actin dynamics: Generating randomness by formin(g) and moving. *Bioarchitecture*. **1**, 165-168, doi:10.4161/bioa.1.4.17314 (2011).
8. Yu, J.H., Crevenna, A.H., Bettenbuhl, M., Freisinger, T., & Wedlich-Soldner, R. Cortical actin dynamics driven by formins and myosin V. *J. Cell. Sci.* **124**, 1533-1541, doi:10.1242/jcs.079038 (2011).
9. Berchtold, D. & Walther, T.C. TORC2 plasma membrane localization is essential for cell viability and restricted to a distinct domain. *Mol. Biol. Cell*. **20**, 1565-1575, doi:10.1091/mbc.E08-10-1001 (2009).
10. Vizcay-Barrena, G., Webb, S.E., Martin-Fernandez, M.L., & Wilson, Z.A. Subcellular and single-molecule imaging of plant fluorescent proteins using total internal reflection fluorescence microscopy (TIRFM). *J. Exp. Bot.* **62**, 5419-5428, doi:10.1093/jxb/err212 (2011).