

Materials List for:

Murine Pancreatic Islet Isolation

Gregory L. Szot¹, Pavel Koudria¹, Jeffrey A. Bluestone¹

¹Diabetes Center, University of California, San Francisco - UCSF

Correspondence to: Gregory L. Szot at gszot@diabetes.ucsf.edu

URL: <http://www.jove.com/video/489>

DOI: [doi:10.3791/489](https://doi.org/10.3791/489)

Materials

Name	Type	Company	Catalog Number	Comments
Dithizone	reagent			Dissolve 5mg Dithizone in 5 ml of DMSO in a 50 ml conical tube. Add 45 ml of normal saline and mix gently. Filter with a 60 cc syringe and a 0.22 "m syringe filter into a new 50 ml tube. Store at 4 OC, until use. Add 1 ml to each 35 mm counting Grid-plate, when needed.
PBS	Buffer			
HBSS	Buffer			
Isolation soln.	Buffer			HBSS containing 10 mM of HEPES, 1 mM of MgCl ₂ , 5 mM of Glucose, pH 7.4
Wash soln.	Buffer			Isolation buffer containing 1mM CaCl ₂
Collagenase P		Boehringer Ingeheim	1213873	
Dissection tools				Hemostat clamp (n=2), Forceps (n=2), Scissors (for skin incision an for taking out pancreas from adjacent tissue)
Syringe				sterile, 5ml, 1/mouse.
30 G needles				sterile, 1/syringe
siliconized glass vials w/ Teflon cap		Fisher Scientific	213-018-54	sterile, 1/mouse
large sink strainer				sterile
1000 ml Beaker				sterile
Gauze pads				4x4, 2/mouse
plastic eye dropper				sterile
Mice	Animal			