

Materials List for:

Bio-layer Interferometry for Measuring Kinetics of Protein-protein Interactions and Allosteric Ligand Effects

Naman B. Shah¹, Thomas M. Duncan¹

¹Department of Biochemistry and Molecular Biology, SUNY Upstate Medical University

Correspondence to: Thomas M. Duncan at duncant@upstate.edu

URL: <https://www.jove.com/video/51383>

DOI: [doi:10.3791/51383](https://doi.org/10.3791/51383)

Materials

Name	Company	Catalog Number	Comments
Octet-RED96	Pall/FortéBio	30-5048	
Bovine Serum Albumin	Sigma	A6003-10G	Fatty Acid free
Biosensor/Streptavidin	Pall/FortéBio	18-5019	Tray of 96 sensors
Microtiter plate	Greiner Bio-one	655209	Black, Polypropylene
Data Acquisition software	Pall/FortéBio	Version 6.4	Newer versions available
Data Analysis software	Pall/FortéBio	Version 6.4	Newer versions available