

Materials List for:

Capture Compound Mass Spectrometry - A Powerful Tool to Identify Novel c-di-GMP Effector Proteins

Benoît-Joseph Laventie¹, Jutta Nesper¹, Erik Ahné², Timo Glatter², Alexander Schmidt², Urs Jenal¹

¹Focal Area Infection Biology, Biozentrum of the University of Basel

²Proteomics Core Facility, Biozentrum of the University of Basel

Correspondence to: Urs Jenal at urs.jenal@unibas.ch

URL: <https://www.jove.com/video/51404>

DOI: [doi:10.3791/51404](https://doi.org/10.3791/51404)

Materials

Name	Company	Catalog Number	Comments
caproBox	caprotec bioanalytics	1-5010-001 (220 V)	UV lamps coupled to a cooling 96-plate cooling block, for the photoactivation
caproMag	caprotec bioanalytics	included in the CCMS Starter Kit	For easy handling of magnetic particles without pipetting
c-di-GMP caproKit	caprotec bioanalytics	upon request	The kit contains the c-di-GMP-capture compound, c-di-GMP (for the competition control), streptavidin coated magnetic beads, capture buffer, and washing buffer
Disposable PD-10 Desalting Columns	GE Healthcare	17-0851-01	
12-tube PCR strips	Thermo Scientific	AB-1114	
UIS250v sonicator with VialTweeter	Hielscher ultrasound technology	UIS250v and VialTweeter	
Miniature French Pressure Cell	Thermo Electron Corporation	FA-003	