

Materials List for:

qPCR Is a Sensitive and Rapid Method for Detection of Cytomegaloviral DNA in Formalin-fixed, Paraffin-embedded Biopsy Tissue

Morgan H. McCoy¹, Kristin Post², Joyashree D. Sen², Hsim Y. Chang², Zijin Zhao¹, Rong Fan¹, Shaoxiong Chen¹, Diane Leland¹, Liang Cheng¹, Jingmei Lin¹

¹Department of Pathology and Laboratory Medicine, Indiana University School of Medicine

²Department of Pathology and Laboratory Medicine, Indiana University Health

Correspondence to: Morgan H. McCoy at momccoy@iupui.edu

URL: <https://www.jove.com/video/51570>

DOI: [doi:10.3791/51570](https://doi.org/10.3791/51570)

Materials

Name	Company	Catalog Number	Comments
1.7 ml Microfuge tubes	Costar	3620	
Serrated forceps	Electron Microscopy Services	62086-1S	Micro Forceps MF-1
Tabletop centrifuges (microfuge)	Eppendorf	5424	
Xylene	Fisher	Fisher Scientific: X3P	1 gallon
Ethanol	Fisher	Fisher Scientific: ET108	96-100%
Microtome	Leica	RM2255	
QIAamp DNA FFPE Tissue Kit	Qiagen	56404	
artus CMV LightCycler PCR ASR	Qiagen	4500025	
CMV LC PCR Supplement Kit	Qiagen	1031873	
LightCycler centrifuge	Roche	75005087	
LightCycler Cooling Block	Roche	10800058001	
LightCycler Capping Tool	Roche	3357317	
Color Compensation Kit	Roche	2158850	
LightCycler 20 µl Capillaries	Roche	1 909 339	
Blades	Sakura	Fisher Scientific: 4689	Accu-Edge low profile