

Materials List for:

Bladder Smooth Muscle Strip Contractility as a Method to Evaluate Lower Urinary Tract Pharmacology

F. Aura Kullmann¹, Stephanie L. Daugherty², William C. de Groat², Lori A. Birder^{1,2}

¹Department of Medicine, Renal division, University of Pittsburgh School of Medicine

²Department of Pharmacology and Chemical Biology, University of Pittsburgh School of Medicine

Correspondence to: F. Aura Kullmann at aura.kullmann@gmail.com

URL: <https://www.jove.com/video/51807>

DOI: [doi:10.3791/51807](https://doi.org/10.3791/51807)

Materials

Name	Company	Catalog Number	Comments
Equipment			
Tissue Bath System with Reservoir	Radnoti, LLC	159920	isolated tissue baths
Warm water recirculator pump	Kent Scientific Corporation	TPZ-749	to keep tissue baths to 37 °C
Computer			
Data Acquisition System	DataQ Instruments	DI-710-UH	To view, record and analyze data
Transbridge Transducer Amplifier	World Precision Instruments	SYS-TBM4M	Transducer amplifier
Grass stimulator	Grass Technologies	Model S88	Stimulator
Anesthesia System	Kent Scientific Corporation	ACV-1205S	To anesthetize the animal
Anesthetizing Box	Harvard Apparatus	500116	To anesthetize the animal
Anesthesia Masks	Kent Scientific Corporation	AC-09508	To anesthetize the animal
Materials and Surgical Instruments			
Sylgard	Dow Corning Corp	184 SIL ELAST KIT	To pin, dissect, & cut tissue
Petri Dish	Corning	3160-152	To dissect/cut tissue
Insect Pins	ENTOMORAVIA Austerlitz Insect Pins	Size 5	To pin tissue
Bench Pad	VWR International	56617-014	Absorbent bench underpads
Rat surgical Kit	Kent Scientific Corporation	INSRATKIT	To remove and dissect tissue
2 Dumont #3 Forceps	Kent Scientific Corporation	INS500064	To remove and dissect tissue
Tissue Forceps	Kent Scientific Corporation	INS500092	To remove and dissect tissue
Scalpel	Kent Scientific Corporation	INS500236	To remove and dissect tissue
Scalpel blade	Kent Scientific Corporation	INS500239	To remove and dissect tissue
Professional Clipper	Braintree Scientific, Inc.	CLP-223 45	To remove fur
Suture Thread	Fine Science Tools	18020-50	Tie tissue
Tissue Clips	Radnoti, LLC	158802	Attach tissue to rod/transducer
1 g weight	Mettler Toledo	11119525	For transducer calibration
Chemicals			
Krebs Solution: Sodium chloride Potassium chloride Monobasic potassium phosphate Magnesium sulfate Dextrose Sodium bicarbonate Calcium chloride Magnesium chloride	Sigma Fisher Fisher Fisher Fisher Sigma EMD Baker	S7653 P217-500 P285-3 M65-500 D16-500 S5761 CX0130-2 2444	To prepare Krebs solution
Isoflurane	Henry Schein	029405	To anesthetize the animal

Oxygen tank	Matheson Tri Gas	ox251	To use with anesthesia system
Carbogen Tank (95% Oxygen; 5% Carbon Dioxide)	Matheson Tri Gas	Moxn00hn36D	To aerate Krebs solutions