

Materials List for:

Clinical Assessment of Spatiotemporal Gait Parameters in Patients and Older Adults

Julia F. Item-Glatthorn¹, Nicola A. Maffiuletti¹

¹Neuromuscular Research Laboratory, Schulthess Clinic

Correspondence to: Julia F. Item-Glatthorn at julia.item@kws.ch

URL: <https://www.jove.com/video/51878>

DOI: [doi:10.3791/51878](https://doi.org/10.3791/51878)

Materials

Name	Company	Catalog Number	Comments
Optogait system (10 meters)	Microgate, Bolzano, Italy		www.optogait.com
Optogait software			www.optogait.com/Support/Downloads
Laptop			
The Optogait system contains the following equipment:			
10 Light-transmitting (T) bars (1 as a first meter)			
10 Light-receiving (R) bars (1 as a first meter)			
18 Caps to connect the bars within a set (9 for T and 9 for R bars) and 2 special caps for the last T and R bar			
1 Camera with its tripod			
1 Cable for connecting the Optogait to the laptop			
1 Cable for connecting the camera to the laptop			
2 Power supplies (one for each set of bars)			