

Materials List for:

Collection, Isolation, and Flow Cytometric Analysis of Human Endocervical Samples

Jennifer A. Juno¹, Genevieve Boily-Larouche¹, Julie Lajoie¹, Keith R. Fowke^{1,2}

¹Department of Medical Microbiology, University of Manitoba

²Department of Community Health Sciences, University of Manitoba

Correspondence to: Julie Lajoie at julie.lajoie1@gmail.com

URL: <https://www.jove.com/video/51906>

DOI: [doi:10.3791/51906](https://doi.org/10.3791/51906)

Materials

Name	Company	Catalog Number	Comments
100 µm cell strainer for 50 ml Falcon tube	BD	352360	CMC processing
RPML 1640	Hyclone	SH30027.01	CMC processing
Fetal bovine serum	Life technology	16000044	CMC processing
Fungizone	Life technology	15290-018	CMC processing
Penicillin/streptomycin	Sigma	P4333-20ml	CMC processing
50 ml Falcon tube	Fisher	14-959-49A	CMC processing
Blood Bank disposable transfer pipette	Fisher	13-711-6M	CMC processing
Cytobrush plus	Cooper surgical	C0121	CMC sampling
Disposable cervical scraper	Quick medical	2183	CMC sampling
15 ml Falcon tube	Fisher	14-959-70c	CVL processing
1.5 ml tube Eppendorf	Fisher	05-402-18	CVL storage
LIVE/DEAD Fixable Cell Stain Kit	Invitrogen	Various	Flow cytometry reagent
Fixation buffer (4% PFA)	BD	554655	Flow cytometry reagent
IgG mouse	Sigma	I8765	Flow cytometry reagent