

Materials List for:

Methods for Skin Wounding and Assays for Wound Responses in *C. elegans*

Suhong Xu¹, Andrew D. Chisholm¹

¹Division of Biological Sciences, Section of Cell and Developmental Biology, University of California, San Diego

Correspondence to: Andrew D. Chisholm at chisholm@ucsd.edu

URL: <https://www.jove.com/video/51959>

DOI: [doi:10.3791/51959](https://doi.org/10.3791/51959)

Materials

Name	Company	Catalog Number	Comments
Agarose	Denville Scientific	CA3510-8	
Plastic Petri plate	Tritech Research	T3308	60 mm
Levamisole	Sigma	L9756	12 mM
Borosilicate glass capillary	World Precision Instruments	1B100F-4	
Needle puller	Sutter Instrument	P-2000	
Worm pick (platinum wire)	Tritech Research	PT-9010	
M9 solution	Home-made		3 g KH ₂ PO ₄ , 6 g Na ₂ HPO ₄ , 5 g NaCl, 1 ml 1 M MgSO ₄ , H ₂ O to 1 liter. Sterilize by autoclaving
Trizol	Invitrogen	15596026	Use 500 ml for 50 worms
iQ SYBR Green	Bio-Rad	1708882	
SuperScript III	Invitrogen	18080-051	
PCR machine	Bio-Rad	C1000/CFX96	
96-well PCR tubes	Bio-Rad	MLL9601	
Image analysis software	Molecular Devices	Metamorph	