

Materials List for:

Measurement of Neurophysiological Signals of Ignoring and Attending Processes in Attention Control

Agatha Lenartowicz¹, Gregory V. Simpson², Samantha R. O'Connell¹, Mark S. Cohen³

¹Department of Psychiatry, University of California Los Angeles

²Attention Research Institute

³Departments of Psychiatry, Radiology, Neurology, Biomedical Physics, Psychology and Bioengineering, University of California Los Angeles

Correspondence to: Agatha Lenartowicz at alenarto@ucla.edu

URL: <https://www.jove.com/video/52958>

DOI: [doi:10.3791/52958](https://doi.org/10.3791/52958)

Materials

Name	Company	Catalog Number	Comments
NetStation Software	Electrical Geodesic, Inc.	version 4.5.1	Alternate recording software may be used.
Matlab Software	The MathWorks, Inc.	7.10.0 (R2010a)	Alternate analysis and presentation software may be used.
PsychToolbox Software	http://psychtoolbox.org/	v3.0.8 (2010-03-06)	Open-source software. Alternate stimulus presentation software may be used.
Netstation Amplifier	Electrical Geodesic, Inc.	300	Alternate amplifier may be used.
EEG Net	Electrical Geodesic, Inc.	HCGSN130	Alternate EEG cap may be used.
Saline-Based Electrolyte (Potassium Chloride)	Electrical Geodesic, Inc.	n/a	Electrolyte used in soaking of net for this high-impedance EEG system. Alternate electrolyte mediate can be used.