

Materials List for:

Use of Galvanic Skin Responses, Salivary Biomarkers, and Self-reports to Assess Undergraduate Student Performance During a Laboratory Exam Activity

Idalis Villanueva¹, Maria Valladares¹, Wade Goodridge¹

¹Department of Engineering Education, Utah State University

Correspondence to: Idalis Villanueva at idalis.villanueva@usu.edu

URL: <https://www.jove.com/video/53255>

DOI: [doi:10.3791/53255](https://doi.org/10.3791/53255)

Materials

Name	Company	Catalog Number	Comments
Salivary Kit	ZRT Laboratories	Female/Male Saliva Profile I	An account needs to be created with the company before sending samples; diurnal samples are recommended for this protocol
Galvanic Skin Response Wrist Sensor and Band	Empatica, Inc.	E3P.v4	M/L fits for the majority of the population although you may need to purchase a smaller size for thin wrists
E-Prime Software	Psychology Software Tools	E-Prime 2.0 Professional	Consider time stamping separately from software