

**Materials List for:**

# **Methods for Acute and Subacute Murine Hindlimb Ischemia**

Michael E. Padgett<sup>1</sup>, Timothy J. McCord<sup>1</sup>, Joseph M. McClung<sup>1</sup>, Christopher D. Kontos<sup>1</sup>

<sup>1</sup>Division of Cardiology, Department of Medicine, Duke University Medical Center

Correspondence to: Christopher D. Kontos at [chris.kontos@duke.edu](mailto:chris.kontos@duke.edu)

URL: <https://www.jove.com/video/54166>

DOI: [doi:10.3791/54166](https://doi.org/10.3791/54166)

## **Materials**

| <b>Name</b>                  | <b>Company</b>          | <b>Catalog Number</b> | <b>Comments</b> |
|------------------------------|-------------------------|-----------------------|-----------------|
| Dumont #5/45 Forceps         | Fine Science Tools      | 11251-35              | Dumoxel         |
| Dumont Style 5 Mini Forceps  | Fine Science Tools      | 11200-14              | Inox            |
| Extra Fine Bonn Scissors     | Fine Science Tools      | 14084-08              |                 |
| 7 - 0 Silk Suture            | Sharp                   | DA-2527N              |                 |
| 5 - 0 Coated Vicryl Suture   | Ethicon                 | J463G                 |                 |
| Graefe Forceps               | Fine Science Tools      | 11053-10              |                 |
| Vannas Spring Scissors       | Fine Science Tools      | 15000-03              |                 |
| Artificial Tears Ointment    | Rugby Laboratories      | 0536-6550-91          |                 |
| Surgical Tape                | 3M                      | 1530-0                |                 |
| Fine Cotton Swabs            | Contec                  | SC-4                  |                 |
| Temperature Controller       | Physitemp               | TCAT-2DF              |                 |
| Ameroid Constrictors         | Research Instruments SW | MMC-0.25 x 1.00-SS    |                 |
| Hot Bead Sterilizer          |                         |                       |                 |
| Deltaphase Isothermal Pad    | Braintree Scientific    | 39DP                  |                 |
| Needle Driver                | Fine Science Tools      |                       |                 |
| Phosphate Buffered Saline    | Gibco                   | 10010-023             |                 |
| Moor LDPI                    | Moor Instruments        | moorLDI2              |                 |
| moorLDI Measurement software | Moor Instruments        | v. 6.0                |                 |
| Hair Removal Cream           | Nair                    |                       |                 |