

Materials List for:

Potentiodynamic Corrosion Testing

Selin Munir¹, Matthew H. Pelletier¹, William R. Walsh¹

¹Surgical and Orthopaedics Research Laboratories, Prince of Wales Clinical School

Correspondence to: Selin Munir at selin.munir@gmail.com

URL: <https://www.jove.com/video/54351>

DOI: [doi:10.3791/54351](https://doi.org/10.3791/54351)

Materials

Name	Company	Catalog Number	Comments
Potentiostat	Metrohm	PGSTAT101	
Ag/AgCl reference electrode, shielded	Metrohm	6.0729.100	
Electrode shaft	Metrohm	6.1241.060	
Polisher Forcipol 1v	Metkon	3602	
Cylindrical flask 700 ml	SciLabware	FR700F	
Reaction lid	SciLabware	MAF2/41	
Dichloromethane	Sigma-Aldrich	MKBR7629V	Use under a fumehood. Wear protective clothing.
Thermo / HAAKE D Series Immersion Circulators	Haake		