

Materials List for:

Removal of *Drosophila* Muscle Tissue from Larval Fillets for Immunofluorescence Analysis of Sensory Neurons and Epidermal Cells

Conrad M. Tenenbaum¹, Elizabeth R. Gavis¹

¹Department of Molecular Biology, Princeton University

Correspondence to: Elizabeth R. Gavis at gavis@princeton.edu

URL: <https://www.jove.com/video/54670>

DOI: [doi:10.3791/54670](https://doi.org/10.3791/54670)

Materials

Name	Company	Catalog Number	Comments
Dumont #5 tweezers	Electron Microscopy Sciences	72701-D	
Micro Scissors, 8 cm, straight, 5 mm blades, 0.1 mm tips	World Precision Instruments	14003	
Sylgard 184 silicone elastomer kit	Dow Corning	3097358-1004	for dissecting plates
Austerlitz insect pins, 0.1 mm	Fine Science Tools	26002-10	
Fostec 8375 light source	Artisan Technology Group	62792-4	
Zeiss Stemi 2000	Carl Zeiss Microscopy		
Vectashield antifade mounting medium	Vector Laboratories	H-1000	for confocal microscopy
Prolong Diamond antifade mountant	Life Technologies	P36970	for structured illumination microscopy
Micro cover glass, 22 x 22 mm, No. 1.5	VWR	48366-227	
Superfrost Plus microscope slides, 25 x 75 x 1.0 mm	Fisherbrand	12-550-15	
Mouse anti-Coracle antibody	Developmental Studies Hybridoma Bank	C615.16	supernatant, dilute 1:50
Mouse anti-Discs large antibody	Developmental Studies Hybridoma Bank	4F3	supernatant, dilute 1:50
Rabbit anti-dsRed antibody	Clontech	632496	dilute 1:1,000
Goat anti-rabbit antibody, Alexa Fluor 568 conjugated	ThermoFisher Scientific	A-11011	dilute 1:1,000
Goat anti-mouse antibody, Alexa Fluor 488 conjugated	ThermoFisher Scientific	A-11001	dilute 1:500