

Materials List for:

Assessing the Effects of Music Listening on Psychobiological Stress in Daily Life

Alexandra Linnemann¹, Jana Strahler¹, Urs M. Nater¹

¹Department of Psychology, Division of Clinical Biopsychology, University of Marburg

Correspondence to: Urs M. Nater at nater@uni-marburg.de

URL: <https://www.jove.com/video/54920>

DOI: [doi:10.3791/54920](https://doi.org/10.3791/54920)

Materials

Name	Company	Catalog Number	Comments
SaliCap Set	IBL International GmbH, a Tecan Group company	RE69985	Sampling tubes for collection of saliva samples to be used in the IBL Saliva Immunoassays
Cortisol Saliva ELISA	IBL International GmbH, a Tecan Group company	RE52611	Enzyme Immunoassay for the quantitative determination of free cortisol in human saliva
Calibrator f.a.s. w/o diluent	Roche Diagnostics	10759350190	Ready-for-use calibrator consisting of a buffered description aqueous solution with biological materials added as required to obtain desired component levels.
Precinorm U	Roche Diagnostics	10171743122	Ready-for-use control
Precipath U	Roche Diagnostics	10171760122	Ready-for-use control
AMY EPS HIT 917 liquid	Roche Diagnostics	11876473 316	R1: α -glucosidase; R2: 4,6-ethylidene-(G7) p-nitrophenyl-(G1)- α ,D-maltoheptaoside
Further materials include typical laboratory utensils, e.g., micropipettes, orbital shaker, vortex mixer, 8-channel micropipettor, wash bottle, automated or semi-automated microtiter plate washing system, precision scale, microtiter plate reader capable of reading absorbance.			
Apple iPod touch, 8 GB, 5 th Generation	Apple Inc.	n/a	mobile diary device
iDialogPad	Mutz Elektronik Entwicklung	n/a	software for programming items occurring in the ambulatory assessment