

Materials List for:

A High-performance Compact Photoacoustic Tomography System for *In Vivo* Small-animal Brain Imaging

Paul Kumar Upputuri¹, Vijitha Periyasamy¹, Sandeep Kumar Kalva¹, Manojit Pramanik¹

¹School of Chemical and Biomedical Engineering, Nanyang Technological University

Correspondence to: Manojit Pramanik at manojit@ntu.edu.sg

URL: https://www.jove.com/video/55811

DOI: doi:10.3791/55811

Materials

Name	Company	Catalog Number	Comments
Pulsed laser diode	Quantel, France	QD-Q1910-SA-TEC	It is the excitation laser source with specifications 803 nm, 1.4 mJ per pulse, 136 ns pulse, 7 kHz maximum, dimentions: 11.0 x 6.0 x 3.6 cm, weight: ~150 gm
Stepper motor with gearbox	LIN Engineering (Servo Dynamics)	Motor: CO-5718-01P, Gearbox: DPL64/1, I = 10 for NEMA 23; power supply PW100-48	To move the detector holder in a circular geometry. Torque: 2.08 N-m, Rotor inertia: 2.6 kg-cm2
Ultrasonic pulser/receiver	Olympus	5072PR	To receive, filter and ampligy the PA signal from UST. Its bandwidth is 35 MHz, and gain is ±59 dB.
Ultrasound Transducer	Olympus	V306-SU-NK-CF1.9IN/Q4200069	Ultrasonic sensors used for photoacoustic detection. Central freqency 2.25 MHz, 0.5 in, Cylindrical focus 1.9 inch
PCIe DAQ (Data acquisition) Card	GaGe	CSE4227/ A6000610/B0E00610	12 bit, 100 Ms/s, 2 channels, 1 Gs on board memory, PCle x16 interface
Rats	In Vivos Pte Ltd, Singapore	NTac:SD, Sprague Dawley / SD	Female, weight 100 ±10g
Acrylic water tank	NTU workshop	Custom-made	It contains the water that acts as an acoustic coupling medium between brain and detector
Circular Scanner	NTU workshop	Custom-made	Scanner is made out of Alluminum
Anesthetic Machine	medical plus pte ltd	Non-Rebreathing Anaesthesia machine with oxygen concentrator.	Supplies oxygen and isoflurane to animal
Pulse Oxymeter portable	Medtronic	PM10N with veterinary sensor	Monitors the pulse oxymetry of the animal
Ultrasound gel	Progress/parker acquasonic gel	PA-GEL-CLEA-5000	Clear ultrasound gel
Data acqusison software	National Instruments Corporation,Austin,TX,USA)	NI LabVIEW 2015 SP1	LabVIEW based program was developed in our laboratory for controlling the stepper motor and acquring the PA singnals from the detector
Data processing software	Matlab (Mathworks, Natick, MA, USA)	Matlab R2012b	Matlab code for reconstruction of PA images was developed in our lab
Temperature controller	LaridTech, MO,USA	MTTC1410	It will constantly control temperature of the PLD
12 V power supply	Voltcraft	PPS-11810	To supply operating voltage for PLD
Variable power supply	BASETech	BT-153	To change the laser output power
Funtion generator	Funktionsgenerator	FG250D	To change the repetetion rate of the PLD. It will provide TTL signal



			to synchronize the DAQ with the laser excitation.
Animal distributor	In Vivos Pte Ltd, Singapore		Animal distributor that supplies small animals for research purpose.
Animal holder	NTU workshop	Custom-made	Used for holding the animal on its abdomen
Breathing mask	NTU workshop	Custom-made	Used along with animal holder to supply anesthesia mixture to the animal
Pentobarbital sodium	Valabarb		Used for euthanizing the animal after the expeirment.
Optical diffuser	Thorlabs	DG10-1500	Used to to make the laser beam homogeneous