

Materials List for:

Central and Divided Visual Field Presentation of Emotional Images to Measure Hemispheric Differences in Motivated Attention

Aminda J O'Hare¹, Ruth Ann Atchley², Keith M Young³

¹Department of Psychology, University of Massachusetts Dartmouth

²Department of Psychology, University of Kansas

³Department of Psychology, University of Minnesota Duluth

Correspondence to: Aminda J O'Hare at aohare@umassd.edu

URL: <https://www.jove.com/video/56257>

DOI: [doi:10.3791/56257](https://doi.org/10.3791/56257)

Materials

Name	Company	Catalog Number	Comments
64-channel Ag-AgCl active electrodes	Cortech Solutions	DA-AT-ESP32102064A/DA-AT-ESP32102064B	EEG electrodes for data collection
ActiveTwo Base System	Cortech Solutions	DA-AT-BCBS	Digitizes and amplifies EEG data at 500 Hz
E-Prime Professional 2.0	Psychology Software Tools	NA	Stimulus presentation software, available at https://www.pstnet.com/eprime.cfm
CURRY 7.0	Compumedics Neuroscan	NA	EEG/ERP data processing and analysis, available at http://compumedicsneuroscan.com/products/by-name/curry/