

Materials List for:

Endotracheal Intubation of Rabbits Using a Polypropylene Guide Catheter.

Krista L Thompson¹, Thomas R Meier¹, Jodi A Scholz¹

¹Department of Comparative Medicine, Mayo Clinic

Correspondence to: Jodi A Scholz at scholz.jodi@mayo.edu

URL: <https://www.jove.com/video/56369>

DOI: [doi:10.3791/56369](https://doi.org/10.3791/56369)

Materials

Name	Company	Catalog Number	Comments
Ketalar Ketamine HCl Inj.	JHP Pharmaceuticals	NDC 42023-113-10	35mg/kg IM in epaxial muscles or quadriceps muscle (preferred)
Anased Xylazine (100mg/ml)	Lloyd Labs	NADA #139-236	5mg/kg IM in epaxial muscles or quadriceps muscle (preferred)
Welch Allyn Laryngoscope Handle	Welch Allyn	60300	Requires 2 size C alkaline batteries
Size 0 or 1 Miller Laryngoscope Blade	Welch Allyn	68040	
2.5 – 4.0 mm ID uncuffed endotracheal tube	Jorgenson	J0149V	Suitable for 2.5-5 kg rabbits 3.5 mm most commonly used
Sovereign polypropylene catheter, 5 Fr x 22"	Medtronic	8890703310	Cut to 8-10 inches Mark the cut end with permanent marker to indicate which end not to insert
Lidocaine 2% HCL Inj	VetOne	L-2000-04	
4-Ply Nonsterile, Bulk Packed Versalon Sponges (Gauze pads)	Covidien	KC9024A	
Tape, Twill, Bleached ½" x 72 yds	Encompass Textiles	48210-025	Used for securing endotracheal tube
HR Lubricating Jelly	HR Pharmaceuticals	201	For lubricating endotracheal tube tip
Great Plains Ballistics: BAAM (Beck Airway Airflow Monitor) 22	BoundTree Medical	21410	For confirmation of successful intubation
Ambu Pediatric Co2 Detector	AMBU INC	172713	For confirmation of successful intubation