

Materials List for:

Live Images of GLUT4 Protein Trafficking in Mouse Primary Hypothalamic Neurons Using Deconvolution Microscopy

Chun Austin Changou^{1,2,3}, Reni Ajoy^{4,5}, Szu-Yi Chou^{4,5,6}

¹The Ph.D. Program for Cancer Biology and Drug Discovery, College of Medical Science and Technology, Taipei Medical University

²Integrated Laboratory, Center of Translational Medicine, Taipei Medical University

³Core Facility, Taipei Medical University

⁴The Ph.D. Program for Neural Regenerative Medicine, College of Medical Science and Technology, Taipei Medical University

⁵Graduate Institute of Neural Regenerative Medicine, College of Medical Science and Technology, Taipei Medical University

⁶TMU research center for Neurotrauma and Neuroregeneration, College of Medical Science and Technology, Taipei Medical University

Correspondence to: Szu-Yi Chou at sichou@tmu.edu.tw

URL: <https://www.jove.com/video/56409>

DOI: [doi:10.3791/56409](https://doi.org/10.3791/56409)

Materials

| Name | Company | Catalog Number | Comments |
|--------------------------------------|---|----------------|--|
| DeltaVision deconvolution microscope | Applied Precision Inc./(GE Healthcare Life Science) | | equipped with 60x/1.42 NA oil immersion objective lens; used for taking live cell images |
| SoftWorX application software | Applied Precision Inc. | | used to control microscope and camera |
| VoloCITY software | PerkinElmer | | used to analyze the images |
| Insulin | Actrapid, Denmark | | |
| Liposome | Invitrogen | 11668019 | Lipofectamine® 2000, used for plasmid DNA transfection |
| GFP control plasmid DNA | | | provided by Professor Samuel Cushman |
| GFP-GLUT4 plasmid DNA | | | provided by Professor Samuel Cushman |
| Anti-mouse Alex-488 | Invitrogen | #A-11001 | Secondary anitbody |
| Anti-rabbit IgG -Alex-568 | Invitrogen | #A-11036 | Secondary anitbody |
| CCL5/RANTES | R&D Systems | 478-MR-025 | 10ng/ml |
| Kimwipes | Kimberly-Clark | #FL42572A | 0.17mm thickness |
| 12 mm Microscope coverglass | Deckglaser | #41001112 | used for immunstaining study |
| micro-dissecting scissors | Klappenecker | | Surgical tool |
| curved-tipped forceps | Ideal-tek | | Surgical tool |
| standard straight-tipped forceps | Ideal-tek | | Surgical tool |
| EndoFree Plasmid Maxi Kit | QIAGEN | 12362 | Purify Endotoxin free plasmid DNA |
| 5 ml pipette | Corning costar | 4487 | dissociate brain tissue |