

Video Article

# Single-Cell Quantification of Protein Degradation Rates by Time-Lapse Fluorescence Microscopy in Adherent Cell Culture

Andrea Brigitta Alber<sup>1</sup>, David Michael Suter<sup>1</sup>

<sup>1</sup>UPSUTER, Institute of Bioengineering (IBI), School of Life Sciences, École Polytechnique Fédérale de Lausanne (EPFL)

Correspondence to: David Michael Suter at [david.suter@epfl.ch](mailto:david.suter@epfl.ch)

URL: <https://www.jove.com/video/56604>

DOI: [doi:10.3791/56604](https://doi.org/10.3791/56604)

Keywords: Developmental Biology, Issue 132, Protein degradation, SNAP-tag, single cells, adherent cell culture, mouse embryonic stem cells, fluorescence microscopy, live cell imaging, pulse labeling, fluorescent dye

Date Published: 2/4/2018

Citation: Alber, A.B., Suter, D.M. Single-Cell Quantification of Protein Degradation Rates by Time-Lapse Fluorescence Microscopy in Adherent Cell Culture. *J. Vis. Exp.* (132), e56604, doi:10.3791/56604 (2018).

## Abstract

Proteins are in a dynamic state of synthesis and degradation and their half-lives can be adjusted under various circumstances. However, most commonly used approaches to determine protein half-life are either limited to population averages from lysed cells or require the use of protein synthesis inhibitors. This protocol describes a method to measure protein half-lives in single living adherent cells, using SNAP-tag fusion proteins in combination with fluorescence time-lapse microscopy. Any protein of interest fused to a SNAP-tag can be covalently bound by a fluorescent, cell permeable dye that is coupled to a benzylguanine derivative, and the decay of the labeled protein population can be monitored after washout of the residual dye. Subsequent cell tracking and quantification of the integrated fluorescence intensity over time results in an exponential decay curve for each tracked cell, allowing for determining protein degradation rates in single cells by curve fitting. This method provides an estimate for the heterogeneity of half-lives in a population of cultured cells, which cannot easily be assessed by other methods. The approach presented here is applicable to any type of cultured adherent cells expressing a protein of interest fused to a SNAP-tag. Here we use mouse embryonic stem (ES) cells grown on E-cadherin-coated cell culture plates to illustrate how single cell degradation rates of proteins with a broad range of half-lives can be determined.

## Video Link

The video component of this article can be found at <https://www.jove.com/video/56604/>

## Introduction

It is well known that cellular proteins undergo extensive turnover, with synthesis and degradation rates being specific for each protein and subject to physiological regulation. Traditionally, protein degradation rates have been measured using bulk methods, such as radioactive pulse chase analysis, or involving protein synthesis inhibitors such as cycloheximide<sup>1</sup>. More recently, stable isotope labeling with amino acids in cell culture (SILAC) in combination with mass spectrometry has been established to quantify protein turnover on a global scale<sup>2</sup>. However, these methods are limited by population averaging, and information about cell-to-cell variability is therefore lost. Furthermore, transient changes in protein degradation that are unsynchronized across the cell population cannot be identified.

Alternatively, protein half-lives can also be determined by fluorescence-based approaches, which often have the advantage of providing single-cell resolution. For example, a photoactivatable green fluorescent protein (paGFP) has been used to determine Oct4 half-life in the early mammalian embryo<sup>3</sup>. Another method to monitor protein decay in living cells is the use of a SNAP-tag in combination with fluorescence time-lapse imaging. The SNAP-tag is a mutant version of the DNA repair enzyme O<sup>6</sup>-alkylguanine DNA-alkyltransferase (AGT) that specifically reacts with benzylguanine (BG) derivatives, which can be coupled to molecular probes<sup>4,5,6</sup>. Therefore, any SNAP-tag fusion protein can be irreversibly labeled with a fluorescent, cell permeable dye. Pulse labeling of a protein of interest fused to the SNAP-tag, followed by washout of the residual dye, allows for monitoring the degradation of the labeled protein population and thus for determining protein half-life. SNAP-tags have been successfully used for pulse-chase labeling of proteins and for determining protein half-lives in adherent cell culture and *in vivo*<sup>5,7,8,9</sup>. A large variety of SNAP-tag substrates covering commonly used fluorescent spectra are commercially available, enabling the selection of the optimal dye for each specific application. Thus, SNAP-tags can also be used for multicolor imaging in combination with other fluorescent fusion proteins or dyes. Cell-impermeable dyes are suitable for labeling of membrane-tethered proteins, whereas cell-permeable dyes are applicable for monitoring both intracellular and membrane-bound proteins. Furthermore, some of these probes exhibit almost no basal fluorescence and only start emitting a strong fluorescent signal upon binding to a SNAP-tag<sup>10</sup>.

This protocol describes how to measure the degradation rates of different proteins of interest in single cells using a SNAP-tag. Here we apply this method to mouse embryonic stem (ES) cells cultured on E-cadherin, but it should be possible to use it with any adherent cultured cell type. We show that pulse labeling of SNAP-tag fusion proteins followed by fluorescence time-lapse imaging allows for determining the single cell half-lives of various proteins of interest and provides an estimate for the cell-to-cell variability of half-lives in a population of cultured cells.

## Protocol

Note: In this study, the E14 ES cell line was used. However, this protocol is directly applicable to any other mouse ES cell line expressing a protein of interest fused to a SNAP-tag, either by tagging the endogenous protein or by using overexpression. For the examples shown in the results section, doxycycline-inducible SNAP-tag fusion cell lines were used (SNAP-tag fused to the following proteins: Nanog, Oct4, Srsf11, or to the fluorescent proteins mOrange2 and sfGFP, and put under the control of a doxycycline-inducible promoter. See<sup>11</sup> for further information on the plasmids used for the generation of the doxycycline-inducible SNAP-tag fusion cell lines). The doxycycline-inducible system can be particularly useful, as it allows for tightly controlling the timing and intensity of the expression of the protein of interest. C-terminal positioning of the SNAP-tag is recommended, as changing the N-terminal amino acid sequence is more likely to alter the half-life of the target protein (N-end rule<sup>12</sup>).

## 1. E-cadherin coating and cell seeding

1. Generate an ES cell line expressing a protein of interest fused to the SNAP-tag<sup>4</sup>.
2. Coat 100 mm cell culture dishes with 4 mL of 0.1% gelatin (diluted in phosphate buffered saline (PBS) without  $\text{Ca}^{2+}/\text{Mg}^{2+}$ ) for 1 h. Remove gelatin and let dry for 1 h. Use immediately or store the coated dishes for up to 2 months.
3. Grow the cell line of interest in ES cell culture medium (Glasgow Minimum Essential Medium, supplemented with 10% ES cell-qualified fetal bovine serum, 2 mM sodium pyruvate, 1% non-essential amino acids, 1% penicillin/streptomycin, 2 mM L-glutamine, 100  $\mu\text{M}$  2-mercaptoethanol, leukemia inhibitory factor (LIF), 3  $\mu\text{M}$  CHIR99021 and 0.8  $\mu\text{M}$  PD184352) on gelatin-coated dishes at 37 °C and 5%  $\text{CO}_2$  for 2 days until reaching a confluence of 10 - 20 Mio cells per dish.  
Note: Here, LIF was produced by transient transfection of HEK 293T cells, followed by supernatant collection and filtering. Each batch of LIF was tested for its potential to maintain pluripotency, but concentrations were not determined. However, recombinant LIF is also commercially available and the concentration commonly used for mouse ES cell culture is 1000 units/mL<sup>13</sup>.
4. Coat a 96-well plate suitable for imaging with 30  $\mu\text{L}$  of recombinant mouse E-cadherin Fc chimera protein (E-cad-Fc) per well (5 ng/ $\mu\text{L}$ , diluted in PBS with  $\text{Ca}^{2+}/\text{Mg}^{2+}$ ). Use stock concentrations of 100 ng/ $\mu\text{L}$ , stored at -80 °C). Avoid extensive pipetting, as E-cad-Fc is very fragile. Incubate at 37 °C for 1.5 h.  
Note: Depending on the microscope, other formats might be used.
5. Aspirate the E-cad-Fc, wash once with 100  $\mu\text{L}$  of PBS with  $\text{Ca}^{2+}/\text{Mg}^{2+}$ , and add 100  $\mu\text{L}$  of pre-warmed ES cell culture medium.
6. Wash the cells of interest with 5 mL of PBS without  $\text{Ca}^{2+}/\text{Mg}^{2+}$ . Aspirate and add 2 mL of Trypsin-EDTA 0.25%. Incubate for 4 min at 37 °C. Add 4 mL of ES cell culture medium, resuspend and spin down at 1000 x g for 4 min. Aspirate the supernatant and resuspend the pellet in 2 mL of fresh ES cell culture medium.  
Note: PBS without  $\text{Ca}^{2+}/\text{Mg}^{2+}$  should be used for washing the cells before trypsinization in order to remove  $\text{Ca}^{2+}$  ions, which are required for cell adhesion. In contrast, for the E-cad-Fc dilution and coating (step 1.4), use PBS with  $\text{Ca}^{2+}/\text{Mg}^{2+}$ , since the adhesion of ES cells to E-cad-Fc coated dishes is  $\text{Ca}^{2+}$ -dependent<sup>14</sup>.
7. Dilute the resuspended cells 1:10 in 1 mL of ES cell culture medium and count, using a counting chamber (load 10  $\mu\text{L}$  for a chamber of 0.1 mm depth).
8. Seed 30,000 cells/ $\text{cm}^2$  in the E-cad-Fc coated imaging plate. Fill up to 200  $\mu\text{L}$  per well with ES cell culture medium. For doxycycline-inducible cell lines, induce with an appropriate dose of doxycycline (optimize dose and timing of the induction in advance. The cell lines used here were treated with 500 ng/mL doxycycline 24 h before imaging). Incubate at 37 °C and 5%  $\text{CO}_2$  and proceed with the pulse labeling and imaging 24 h after cell seeding.

## 2. Pulse labeling of the SNAP-tag

Note: For protein decay experiments it is crucial to use an adequate SNAP dye concentration. The concentration should be high enough to yield a bright signal in the beginning of the time-lapse, as the fluorescence will decrease over time. However, using too high dye concentrations might cause residual dye being left in the medium or in the cells even after washing. The free dye might subsequently bind to newly produced SNAP-tag molecules over the course of the movie, which will distort the decay curve. The observed fluorescence signal will depend on the properties of the dye, the cell line used, as well as the expression level of the corresponding protein. Therefore, it is crucial to optimize the dye concentration by testing different dilutions, starting from the dilution suggested for live cell imaging by the manufacturer. For this study, a far-red fluorescent substrate was used. An optimal concentration of 12 nM was determined for doxycycline-inducible overexpression cell lines.

1. Dilute the SNAP dye to the appropriate concentration in pre-warmed ES cell culture medium. Use a pipette to gently remove the medium from the cells previously seeded in the 96-well plate and add 50  $\mu\text{L}$  of diluted SNAP dye per well. Incubate for 30 min at 37 °C.
2. Aspirate the diluted dye with a pipette and wash 3x with 200  $\mu\text{L}$  of pre-warmed PBS (without  $\text{Ca}^{2+}/\text{Mg}^{2+}$ ). Add 200  $\mu\text{L}$  of ES cell culture medium and incubate for 15 min at 37 °C.
3. Repeat the previous washing steps twice more.  
Note: Extensive washing is important in order to remove any unbound dye. The repeated washing steps with PBS remove the residual extracellular dye in the medium, whereas the 15 min incubations ensure robust labeling of the SNAP-tag fusion proteins prior to starting the imaging experiment. However, perform the washing steps by gentle pipetting and do not use an automatic aspirator, as cells seeded on E-cad-Fc tend to detach easily.
4. Add 200  $\mu\text{L}$  of imaging medium. To reduce the fluorescence background, use phenol red free medium (supplemented with 10% ES cell-qualified fetal bovine serum, 2 mM sodium pyruvate, 1% non-essential amino acids, 1% penicillin/streptomycin, 2 mM L-glutamine, 100  $\mu\text{M}$  2-mercaptoethanol, LIF, 3  $\mu\text{M}$  CHIR99021 and 0.8  $\mu\text{M}$  PD184352).

### 3. Time-lapse microscopy

1. Use a microscope suitable for live imaging, allowing controlled temperature and CO<sub>2</sub>. Set the temperature to 37 °C and the CO<sub>2</sub> to 5% prior to use and allow to equilibrate for 1 - 2 h.
2. Place the plate into the microscope and find spots suitable for imaging. Select spots where the cells are evenly distributed, but not too dense, as this will facilitate the analysis. Adjust the number of selected spots to the number of cells required for the analysis. In this study, 3 - 5 spots per cell line were recorded.
3. Select the objective. A 20X objective is appropriate for the size of ES cells.
4. Select the illumination settings for the fluorescent channel(s) to be recorded. Adjust the laser power and exposure according to the intensity of the fluorescence signal from the cells of interest. Make sure to obtain a strong initial signal, since it will decrease over the course of the movie, but avoid laser powers higher than 30% in order to minimize phototoxicity and photobleaching. For this study, an excitation filter of 632/22 nm (wavelength/bandpass), emission filter of 679/34 nm (wavelength/bandpass), laser power of 10% and exposure time of 100 ms were used.
5. Select the imaging intervals and duration of the time-lapse experiment. For proteins with expected half-lives of 2 - 20 h, choose acquisition times of 12 - 24 h with time intervals of 15 min. For shorter-lived proteins, reduce the intervals and acquisition times, for longer-lived proteins increase accordingly.
6. Start imaging.

### 4. Image processing and analysis

1. Collect the acquired images as a stack in tiff format. For further processing and analysis, use the FIJI software<sup>15</sup>. If the time-lapse data is not collected as a stack by the microscope software, open all frames of the time-lapse experiment in the software (either by clicking **File | Open...** or drag and dropping the corresponding files to the toolbar) and click **Image | Stacks | Images to Stack**.
2. Use the subtract background function to remove the background from all pictures in the stack. To do so, click **Process | Subtract Background**. Select the rolling ball radius in the pop-up window. Use a radius that is at least the size of the cells imaged. To apply the background subtraction to all cells in the stack, select **Yes** in the **Process stack?** pop-up window.
3. Select a cell of interest and draw a region of interest (ROI) around it using the toolbar. For nuclear signals, use the **Oval** selection by first clicking on the corresponding tab in the toolbar and then clicking on the nucleus of interest and adjusting the oval around it. For cytoplasmic signals or very dense regions use the freehand selection to be able to closely follow the cell outlines. Avoid including too large regions of background, even though it is not necessary to follow the cell outlines precisely due to the subsequent subtraction of all background pixels (step 4.8 - 4.9).
4. Add the region of interest to the ROI manager by clicking **Analyze | Tools | ROI manager | Add**, or by pressing **T** on the keyboard.
5. Proceed to the next frame by moving the tab in the stack window or by pressing shift + < on the keyboard, then repeat the previous actions. Follow the cell of interest throughout the course of the movie and add the ROI of each frame to the ROI manager. Save the final ROI set for each tracked cell by clicking **More | Save...** in the ROI manager.  
Note: After cell divisions, track both daughter cells separately and sum their integrated intensities after background subtraction (see steps 4.6 - 4.9) to account for protein dilution during cell division. Save the ROI sets for both daughter cells.
6. Once the cell of interest is tracked throughout the movie, click **Measure**, to obtain values for the mean intensity and the area of the ROIs. Copy the obtained values into an electronic spreadsheet program and calculate the raw integrated intensity of the cell for each time-point as follows:  
$$I_C = Mean_C * Area_C$$
where  $Mean_C$  is the mean intensity and  $Area_C$  the area of the corresponding ROI.
7. To estimate the local background, add a ROI close to the cell of interest for each time frame. Avoid including any cellular fluorescence signal. Use a circular ROI that is roughly the size of a cell and move or shrink it accordingly if needed, e.g. if neighboring cells interfere. Proceed as described previously with the cellular ROI set in order to obtain a background ROI set and copy the measured intensity values to the spreadsheet.
8. To obtain a background-corrected value for the integrated intensity of the cell, first calculate the integrated intensity of the background for each time-point:  
$$I_{BG} = Mean_{BG} * Area_C$$
where  $Mean_{BG}$  is the mean intensity of the background signal and  $Area_C$  is the area of the ROI encircling the cell. Do not use the area of the background ROI, unless both areas have the same size.
9. Calculate the final background-subtracted integrated intensity of the cell for each time point:  
$$I_{C-BG} = I_C - I_{BG}$$
10. To normalize the single cell decay curves, divide the intensity value of each time-point by the intensity value of the first time-point.  
Note: The curve fitting (see step 4.11) can either be performed on each single cell or on the population average. A normalization is required if a population based average is calculated in order to avoid biases from different fluorescence intensities between cells. The normalization thus ensures that each cell contributes with the same weight to the final decay curve. In addition, the normalization can be useful to visualize the single cell decays independently of their absolute fluorescence intensities (see **Figures 3A-3C**). However, if only the single-cell half-lives are determined and the data is not averaged, the normalization step may be omitted and the curve fitting can be directly performed on the raw data.
11. For the estimation of the protein half-life, use a curve fitting tool. In this study, the MATLAB curve fitting toolbox 3.4.1 was used, which is located by default in the **Apps** section of the MATLAB user interface. Import the fluorescence intensity and time values from the electronic spreadsheet into MATLAB by clicking on the **Import Data** tab. Open the curve fitting toolbox and select the time points and fluorescence decay data in the **X data** and **Y data** tab. Choose **Custom Equation** in the curve fitting tab and enter the equation for an exponential decay:  
$$f(t) = a * e^{-bt}$$
where  $f(t)$  is the fluorescence intensity at a given time-point,  $a$  the initial intensity and  $b$  the decay rate. In the **Fit Options...** tab, select 0 for the lower limit of both  $a$  and  $b$ . The estimated values for  $a$  and  $b$  will then appear in the results window. Calculate the half-life as follows:

$$t_{1/2} = \frac{\ln(2)}{b}$$


## Representative Results

The described protocol provides an estimate of the cell-to-cell variability in half-life for any given protein fused to a SNAP-tag. The use of recombinant E-cadherin-Fc for coating of the imaging plate allows for single cell resolution in ES cells, which otherwise grow in colonies. Single cells can be tracked separately throughout the course of the movie ( **Figure 1A**).

In order to determine the protein half-life for each single cell, it is crucial to measure the integrated, background-subtracted SNAP-tag fluorescence signal over time ( **Figure 1B**), with summing up the integrated intensities of both daughter cells in case of divisions. This results in an exponential decay curve for each cell, from which the decay rate and thus the half-life can be extracted by curve fitting ( **Figure 1C**). Importantly, if an average decay curve is calculated, the single cell traces should be normalized to the first frame to ensure that each cell has the same weight on the average, despite putative differences in initial fluorescence intensity between cells.

The adequate dye concentration depends on the type of dye and cell line used, and should thus be optimized for each experiment. To illustrate this, **Figure 2** shows the decay curves of different concentrations of the far-red fluorescent substrate used for all experiments shown here, tested on a doxycycline inducible SNAP-tag fusion cell line. An initial dye concentration of 3  $\mu\text{M}$  was chosen according to the manufacturer's instructions and a serial dilution of 1:3 was performed until reaching a concentration of 1.4 nM. For the two highest concentrations the signal does not decrease, whereas the observed decay is exponential for concentrations below 111 nM. An optimal concentration of 12 nM of dye was chosen for further experiments, since it ensured minimal amounts of residual dye left in the medium after washing, but the signal was still bright enough to observe clear decay curves for a variety of doxycycline inducible cell lines.

Figures 3A-3C show representative results, including the single cell decays and population averages, for 3 proteins with different average half-lives: the pluripotency associated transcription factors Nanog (average half-life 2.9 h) and Oct4 (average half-life: 4.8 h), and the much longer lived pre-mRNA processing factor Srsf11 (average half-life 25.3 h). The boxplots ( **Figure 3D**) represent the half-lives obtained from the single cell decay curves by curve fitting and illustrate the heterogeneity of half-lives for the three proteins.




**Figure 1: Time-lapse acquisition and analysis**

A: Series of representative images showing the decay of fluorescence in cells pulse labeled with SNAP dye. Cell line: doxycycline inducible TRE3G-Oct4-SNAPtag, induced with 500 ng/mL doxycycline 7 h before imaging.


B: Movie analysis and calculation of integrated intensity.  $I_C$ : Integrated intensity of the cell,  $I_{BG}$ : Integrated intensity of the background,  $I_{C-BG}$ : background-corrected integrated intensity.

C: Exponential curve fitting performed in MATLAB. Example of a single cell decay trace, the corresponding fit and the resulting estimates for the decay parameters  $a$  and  $b$ . The half-life is calculated as follows:  $t_{1/2} = \ln(2)/b$ . [Please click here to view a larger version of this figure.](#)



**Figure 2: Optimization of the dye concentration.**

Fluorescence decay of cells pulse labeled with different concentrations of SNAP dye. A serial dilution of the dye starting from 3  $\mu$ M was performed. The traces are averages of 5 cells each, normalized to the first frame and tracked for 12 h. Cell line: doxycycline inducible TRE3G-sfGFP-mOrange2-SNAPtag-PEST, induced with 500 ng/mL doxycycline 24 h before imaging. [Please click here to view a larger version of this figure.](#)



**Figure 3: Single cell variability of half-lives for different SNAP-tagged candidate proteins.**

A-C: Single-cell decays (light blue) and population averages (dark blue) for doxycycline inducible Nanog-, Oct4- and Srsf11-SNAP-tag fusion proteins.

D: Half-lives of individual cells, which were calculated by exponential curve fitting for each single cell. The boxplots are displayed in  $\log_2$  scale. The boxes display the interquartile range (IQR), the black line indicates the median. The whiskers represent the minimum and maximum values, excluding the outliers (black dots), which are defined as values deviating more than 1.5x from the IQR. N = 20 cells each. [Please click here to view a larger version of this figure.](#)

## Discussion

The most crucial step when using a SNAP-tag to monitor protein decay is to ensure that no residual unbound dye is left in the medium or in the cells after washing, as otherwise it might bind to newly produced SNAP-tag molecules later in the course of the experiment and thereby compromise the decay curve. This is on one hand achieved by carefully performing all the described washing steps. On the other hand, the dye concentration should be kept as low as possible, while still being in a range that allows for obtaining a good signal to noise ratio. As shown in **Figure 2**, the dye concentration should be optimized for each experiment by testing different dilutions. In our case, the doxycycline inducible cell lines have relatively high expression levels, thus allowing to use a low dye concentration.

This protocol describes a general approach to quantify protein degradation in single cells by time-lapse fluorescence imaging and various aspects of the protocol can be modified, depending on the experimental conditions and the goal of the corresponding experiment. Protein degradation mostly follows a first order exponential decay. For the protein decays shown here, the  $R^2$  values of the fits were fairly high (0.95-0.99) and using multi-exponentials has not significantly improved the fits. However, if the obtained results from single exponential fits are not satisfying, fitting multi-exponentials might be considered, which would imply the existence of multiple labeled protein subpopulations



decaying with different rates. Furthermore, our image analysis pipeline involves substantial manual work, which is only feasible if limited numbers of cells are analyzed. For higher numbers of cells, a more automated approach should be considered. Various automated tools have been developed in the past years, which allow for segmenting and tracking larger numbers of cells<sup>16,17</sup>. However, the movement of the cells and frequent cell divisions are challenging and prone to cause segmentation errors. In addition, when imaging protein decays, the fluorescence signal will decrease over time, making the selection of a threshold for segmentation even more difficult. Thus, tracking tools should be carefully tested before being implemented for this protocol, as segmentation errors can distort the data considerably. In addition, a different background subtraction method may be considered, depending on the quality of the acquired images. Here we use FIJI's rolling ball correction, which is suitable for evenly illuminated images with rather sparsely seeded cells. However, especially if images suffer from uneven illumination, more specialized background subtraction methods might be applied<sup>18,19</sup>.

One limitation to our approach is the duration of the time-lapse recordings that can be analyzed. For relatively short-lived proteins, an acquisition time of about 12 h is sufficient. For longer-lived proteins, as for example Srsf11, performing a curve fitting on a time range of 12 h is possible. However, it is important to keep in mind that the exponential fits might be less precise for longer-lived proteins, for which the acquisition time and the imaging intervals should be increased. If fast cycling cells are being imaged, longer acquisition times will complicate the analysis, due to the increasing number of daughter cells to be analyzed. Again, in this case, a more automated analysis pipeline might be useful. In addition, when using SNAP-tag fusion proteins to determine protein decay rates, the possibility of the SNAP-tag altering the half-life of the endogenous protein should be considered. In the case of the cell lines used here we see no evidence for such effects, as the determined half-lives of Oct4 and Nanog closely match published values<sup>3,20,21,22</sup>. To further rule out an effect of the SNAP-tag on protein degradation, especially when no data on half-life is available, one option would be to block protein synthesis by cycloheximide and to perform western blots at different time points using an antibody recognizing the endogenous protein, and directly compare the decay of the SNAP-tagged candidate protein with the untagged endogenous version.

Most previously published data on protein half-life is based on bulk biochemical experiments, mainly time-course western blot experiments upon inhibition of protein synthesis by cycloheximide. However, this approach does not provide single cell resolution and its temporal resolution is low due to the limited number of time-points that can be taken into account. In contrast, our single-cell approach does not require the use of protein synthesis inhibitors and can identify cell-to-cell variability in protein degradation rates. It also allows for identifying potential subpopulations of cells with different half-lives, which might be particularly interesting for the investigation of heterogeneously expressed proteins. A further advantage of our method is the possibility to *in silico* synchronize the cells and therefore monitor transient changes in half-life, for example throughout the cell cycle. However, the SNAP-tag technology requires the generation of a corresponding SNAP-tagged cell line, which might be its major drawback. Importantly, our approach can also be applied to endogenously tagged proteins<sup>9</sup>. As mentioned above we do not observe any evidence of the SNAP-tag altering the endogenous protein half-life for the proteins we have studied, suggesting that the accuracy of the SNAP-tag approach is in the same range as other methods.

The method presented here will be useful for various applications. Most importantly, it provides the option to study cell-to-cell variability in half-lives for any SNAP-tag fusion protein of interest. There are several other tags available (such as the CLIP-tag<sup>23</sup> or the Halo-tag<sup>24</sup>) whose working principle is the same. Pulse labeling and washing of the residual dye will lead to a decay curve, which allows for determining protein half-life, whereas leaving the dye in the medium (in case of a ligand that is non-fluorescent unless bound to the tag) enables long-term imaging of a protein of interest. Due to the large selection of available fluorescent ligands, the SNAP-tag can thus be combined with other tags or with multiple other fluorescently labeled molecules, allowing for applying the technique in various contexts.

## Disclosures

The authors have nothing to disclose.

## Acknowledgements

Time-lapse microscopy experiments were performed at the Biomolecular Screening Facility (BSF), EPFL. We thank Marc Delachaux (Service Audiovisuel, EPFL) for the videography and editing of the movie.

## References

1. Zhou, P. Determining Protein Half-Lives. *Signal Transduct Protoc.* 67–77 (2004).
2. Schwanhäusser, B. *et al.* Global quantification of mammalian gene expression control. *Nature*. **473** (7347), 337–342 (2011).
3. Plachta, N., Bollenbach, T., Pease, S., Fraser, S. E., & Pantazis, P. Oct4 kinetics predict cell lineage patterning in the early mammalian embryo. *Nature Cell Biol.* **13** (2), 117–123 (2011).
4. Keppler, A., Gendreizig, S., Gronemeyer, T., Pick, H., Vogel, H., & Johnsson, K. A general method for the covalent labeling of fusion proteins with small molecules *in vivo*. *Nature Biotechnol.* **21** (1), 86–89 (2002).
5. Keppler, A., Pick, H., Arrivoli, C., Vogel, H., & Johnsson, K. Labeling of fusion proteins with synthetic fluorophores in live cells. *Proc Natl Acad Sci U S A*. **101** (27), 9955–9959 (2004).
6. Gronemeyer, T., Chidley, C., Juillerat, A., Heinis, C., & Johnsson, K. Directed evolution of O6-alkylguanine-DNA alkyltransferase for applications in protein labeling. *Protein Eng Design Select.* **19** (7), 309–316 (2006).
7. Jansen, L. E. T., Black, B. E., Foltz, D. R., & Cleveland, D. W. Propagation of centromeric chromatin requires exit from mitosis. *J Cell Biol.* **176** (6), 795–805 (2007).
8. Bojkowska, K. *et al.* Measuring In Vivo Protein Half-Life. *Chem Biol.* **18** (6), 805–815 (2011).
9. Mandic, A., Streibinger, D., Regali, C., Phillips, N. E., & Suter, D. M. A novel method for quantitative measurements of gene expression in single living cells. *Methods*. **120**, 65–75 (2017).

10. Komatsu, T. *et al.* Real-Time Measurements of Protein Dynamics Using Fluorescence Activation-Coupled Protein Labeling Method. *J Am Chem Soc.* **133** (17), 6745–6751 (2011).
11. Deluz, C. *et al.* A role for mitotic bookmarking of SOX2 in pluripotency and differentiation. *Genes Dev.* **30** (22), 2538–2550 (2016).
12. Varshavsky, A. The N-end rule pathway and regulation by proteolysis. *Protein Sci.* **20** (8), 1298–1345 (2011).
13. Tamm, C., Pijuan Galitó, S., & Annerén, C. A Comparative Study of Protocols for Mouse Embryonic Stem Cell Culturing. *PLoS ONE.* **8** (12), e81156 (2013).
14. Nagaoka, M. *et al.* E-Cadherin-Coated Plates Maintain Pluripotent ES Cells without Colony Formation. *PLoS ONE.* **1** (1), e15 (2006).
15. Schindelin, J. *et al.* Fiji: an open-source platform for biological-image analysis. *Nature Meth.* **9** (7), 676–682 (2012).
16. Hilsenbeck, O. *et al.* Software tools for single-cell tracking and quantification of cellular and molecular properties. *Nat Biotech.* **34** (7), 703–706 (2016).
17. Blanchoud, S., Nicolas, D., Zoller, B., Tidin, O., & Naef, F. CAST: An automated segmentation and tracking tool for the analysis of transcriptional kinetics from single-cell time-lapse recordings. *Methods.* **85**, 3–11 (2015).
18. Peng, T. *et al.* A BaSiC tool for background and shading correction of optical microscopy images. *Nature Comm.* **8** (2017).
19. Smith, K. *et al.* CIDRE: an illumination-correction method for optical microscopy. *Nature Methods.* **12** (5), 404–406 (2015).
20. Chae, H.-D., Lee, M.-R., & Broxmeyer, H. E. 5-Aminoimidazole-4-carboxamide Ribonucleoside Induces G1/S Arrest and Nanog Downregulation via p53 and Enhances Erythroid Differentiation. *Stem Cells.* **30** (2), 140–149 (2012).
21. Lin, Y. *et al.* Reciprocal Regulation of Akt and Oct4 Promotes the Self-Renewal and Survival of Embryonal Carcinoma Cells. *Mol Cell.* **48** (4), 627–640 (2012).
22. Wei, F., Scholer, H. R., & Atchison, M. L. Sumoylation of Oct4 Enhances Its Stability, DNA Binding, and Transactivation. *J Biol Chem.* **282** (29), 21551–21560 (2007).
23. Gautier, A. *et al.* An Engineered Protein Tag for Multiprotein Labeling in Living Cells. *Chem Biol.* **15** (2), 128–136 (2008).
24. Los, G. V. *et al.* HaloTag: A Novel Protein Labeling Technology for Cell Imaging and Protein Analysis. *ACS Chem Biol.* **3** (6), 373–382 (2008).