

Materials List for:

Luciferase Complementation Imaging Assay in *Nicotiana benthamiana* Leaves for Transiently Determining Protein-protein Interaction Dynamics

Kaiwen Sun^{*1}, Yuyu Zheng^{*1}, Ziqiang Zhu¹

¹Jiangsu Key Laboratory for Biodiversity and Biotechnology, College of Life Sciences, Nanjing Normal University

*These authors contributed equally

Correspondence to: Ziqiang Zhu at zqzhu@icloud.com

URL: <https://www.jove.com/video/56641>

DOI: [doi:10.3791/56641](https://doi.org/10.3791/56641)

Materials

Name	Company	Catalog Number	Comments
Transformation solution			
10 mM Morpholineethanesulfonic acid	VETEC	V900336	
27.8 mM Glucose	VETEC	V900392	
10 mM MgCl ₂ •6H ₂ O	VETEC	V900020	
150 μM Acetosyringone	ALDRICH	D134406	
pH 5.7			
Luciferin working buffer			
5 mM Luciferin potassium salt	GOLD BIOTECHNOLOGY	LUCK-100	
0.025% Triton X-100	VETEC	V900502	
H ₂ O to 10 ml			