

Materials List for:

Gait Analysis of Age-dependent Motor Impairments in Mice with Neurodegeneration

Christine M. Rostosky^{1,2}, Ira Milosevic^{1,2}

¹European Neuroscience Institute (ENI)

²University Medical Center Göttingen (UMG)

Correspondence to: Ira Milosevic at imilose@gwdg.de

URL: <https://www.jove.com/video/57752>

DOI: [doi:10.3791/57752](https://doi.org/10.3791/57752)

Materials

Name	Company	Catalog Number	Comments
DigiGait	Mouse Specifics, Inc., Framingham, Massachusetts, USA		DigiGait Imager and Analysis Software are included with the hardware
non-transparent blanket or dark cloth			cover the test chamber to reduce the animal's feeling of exposure/ stress
balance	e.g. Satorius		balance with 0.1 g accuracy and a maximum load of at least 100 g
red finger paint	e.g. Kreul or Staedtler		for increasing the contrast between paws and animal's body
small paint brush			soft brush to apply finger paint to the animal paws
diluted detergent			for cleaning
disinfectant, e.g. Meliseptol or 70% ethanol	e.g. B.Braun		for disinfection