

Video Article

# Oral Intubation of Adult Zebrafish: A Model for Evaluating Intestinal Uptake of Bioactive Compounds

Jie Ji<sup>1,2</sup>, Rosemary Thwaite<sup>1,2</sup>, Nerea Roher<sup>1,2</sup>

<sup>1</sup>Institute of Biotechnology and Biomedicine (IBB-UAB), Universitat Autònoma de Barcelona

<sup>2</sup>Department of Cell Biology, Animal Physiology and Immunology, Universitat Autònoma de Barcelona

Correspondence to: Nerea Roher at [nerea.roher@uab.cat](mailto:nerea.roher@uab.cat)

URL: <https://www.jove.com/video/58366>

DOI: [doi:10.3791/58366](https://doi.org/10.3791/58366)

Keywords: Immunology and Infection, Issue 139, Zebrafish, oral intubation, intestine, mucosal immunity, cytometry, confocal imaging, qPCR

Date Published: 9/27/2018

Citation: Ji, J., Thwaite, R., Roher, N. Oral Intubation of Adult Zebrafish: A Model for Evaluating Intestinal Uptake of Bioactive Compounds. *J. Vis. Exp.* (139), e58366, doi:10.3791/58366 (2018).

## Abstract

Most pathogens invade organisms through their mucosa. This is particularly true in fish as they are continuously exposed to a microbial-rich water environment. Developing effective methods for oral delivery of immunostimulants or vaccines, which activate the immune system against infectious diseases, is highly desirable. In devising prophylactic tools, good experimental models are needed to test their performance. Here, we show a method for oral intubation of adult zebrafish and a set of procedures to dissect and prepare the intestine for cytometry, confocal microscopy and quantitative polymerase chain reaction (qPCR) analysis. With this protocol, we can precisely administer volumes up to 50  $\mu$ L to fish weighing approximately 1 g simply and quickly, without harming the animals. This method allows us to explore the direct *in vivo* uptake of fluorescently labelled compounds by the intestinal mucosa and the immunomodulatory capacity of such biologics at the local site after intubation. By combining downstream methods such as flow cytometry, histology, qPCR and confocal microscopy of the intestinal tissue, we can understand how immunostimulants or vaccines are able to cross the intestinal mucosal barriers, pass through the lamina propria, and reach the muscle, exerting an effect on the intestinal mucosal immune system. The model could be used to test candidate oral prophylactics and delivery systems or the local effect of any orally administered bioactive compound.

## Video Link

The video component of this article can be found at <https://www.jove.com/video/58366/>

## Introduction

The goal of this article is to describe in depth a straightforward method for oral intubation of zebrafish, along with useful associated downstream procedures. Oral intubation using zebrafish has become a practical model in the study of infectious disease dynamics, oral vaccine/immunostimulant, drug/nanoparticle uptake and efficacy, and intestinal mucosal immunity. For example, zebrafish oral intubation has been used in the study of *Mycobacterium marinum* and *Mycobacterium peregrinum* infection<sup>1</sup>. Lovmo *et al.* also successfully used this model to deliver nanoparticles and *M. marinum* to the gastro-intestinal tract of adult zebrafish<sup>2</sup>. In addition, Chen *et al.* used zebrafish oral intubation to show that drugs encapsulated by nanoparticles, when administered *via* the gastro-intestinal tract, were transported across the blood brain barrier<sup>3</sup>. These authors performed intubation based on the gavage method described by Collymore *et al.*<sup>4</sup> with some modifications. However, they did not provide a highly detailed protocol describing the oral intubation procedure. Here, we present a method for oral intubation of adult zebrafish building on Collymore *et al.*<sup>4</sup> We further include the preparation of the intestine for relevant downstream analysis by cytometry, confocal microscopy and qPCR.

The intestine and particularly its mucosa is the first-line of defense against infection and the primary site of nutrient uptake<sup>5</sup>. When the epithelial cells and antigen-presenting cells within mucosal barriers perceive danger signals, an immediate innate immune response is triggered. Next, the highly specific adaptive immune response is established by T and B lymphocytes<sup>6,7</sup>. Development of oral vaccines is a current focus area in vaccinology. Such vaccines would be an effective tool to protect the organism at exposed sites due to the specific response of immune cells in the mucosa-associated lymphoid tissues (MALT)<sup>8,9</sup>. In aquaculture, mucosal vaccines have obvious advantages compared to injectable vaccines. They are practical for mass vaccination, less labor-intensive, are less stressful to the fish, and can be administered to young fish. Nevertheless, mucosal vaccine candidates must reach the second gut segment without being denatured in the oral environment. They also must cross mucosal barriers in order to gain access to antigen presenting cells (APCs) to induce local and/or systemic responses<sup>10</sup>. Hence, testing of the mucosal uptake achieved by candidate oral antigens and their delivery systems, as well as the immune response evoked, is essential in the development of oral vaccines.

In a biomedical context, developing a model to test biological effects of compounds after oral intubation is of growing interest. Many of the anatomical and physiological features of the intestine are conserved between bilaterian lineages, with mammals and bony fishes<sup>11</sup>. This oral intubation model connected to downstream analysis can be a tool to provide insights into human biology, as well as a testing ground for biologics or other compounds *in vivo*.

The oral intubation protocol can be performed by one operator, e.g., successfully administering up to 50  $\mu$ L of the protein nanoparticle suspension to fish weighing 1 g, with a high survival rate. The procedure is simple to set up and quick; 30 fish can be intubated in 1 h. The protocol for intestine preparation is key to providing quality cell and tissue samples for subsequent analysis. Examples of downstream results are given which show the protocol's usefulness in obtaining data related to intestinal uptake and in isolating quality RNA for qPCR. The protocol would be of great use to those needing a suitable model to test the dynamics of oral prophylactics or other compounds in the intestine.

## Protocol

All experimental procedures involving zebrafish (*Danio rerio*) were authorized by the Ethics Committee of the Universitat Autònoma de Barcelona (CEEH number 1582) in agreement with the International Guiding Principles for Research Involving Animals (EU 2010/63). All experiments with live zebrafish were performed at 26–28 °C.

## 1. Preparing the Equipment for Oral Intubation

1. Place approximately 1 cm of a fine silicone tube on a 31 G Luer lock needle to cover the needle tip.
  2. Cut a 10  $\mu$ L sterile filter pipette tip (around 2 cm), take the finer end and place it over the silicone tube as a sheath. Make sure the pipette extends beyond the tip of the needle to avoid injuring the animal.
  3. Attach the needle to a 100  $\mu$ L Luer lock syringe.
- NOTE: Always rinse with ethanol and then phosphate buffered saline (PBS, see **Materials**) thoroughly between treatments.

## 2. Solutions Required

1. Prepare 150 mg/L (for anesthesia) or 300 mg/L (for euthanasia) of ethyl 3-aminobenzoate methanesulfonate (MS-222) solution with water from the aquarium where the zebrafish are maintained. Fill a small tank with 1 L of anesthetic solution and keep it aerated.
  2. Fill another small tank with 1 L of aquarium water without MS-222 for fish recovery and keep it aerated.
  3. Make 50 mL of 1x PBS from a 10x sterile stock solution.
  4. For cytometry analysis/intestinal cell isolation, prepare enough fresh 0.15% collagenase Type IV solution for 1 mL per fish from a stock solution or from powder in Dulbecco's modified eagle medium (DMEM) with 1% v/v penicillin and streptomycin (see materials). Make aliquots (1 per fish) of 1 mL in 2 mL centrifuge tubes. Keep at 4 °C until 30 min before the dissection step.
  5. For confocal microscopy/sample fixation, prepare 50 mL of fresh 4% paraformaldehyde (PFA) solution in PBS or thaw a stock solution from -20 °C freezer in a fume hood.
- CAUTION: PFA is toxic. Please read the material safety data sheet before working with it. Gloves and safety glasses should be worn, and always leave solutions inside a fume hood.

## 3. Preparing the Fluorescent Nanoparticle Suspension

1. Label the protein nanoparticle with Atto-488 NHS ester (see **Table of Materials**) or an appropriate fluorescent dye according to the manufacturer's instructions.
2. Resuspend the nanoparticles in 0.1 M sodium bicarbonate buffer at the concentration of 2 mg/mL.
3. Dissolve the Atto 488 NHS ester in amine-free dimethyl sulfoxide (DMSO) at 2 mg/mL. Keep an aliquot of 10  $\mu$ L to check labeling efficiency (step 3.7–3.8).
4. Mix the nanoparticles and the Atto 488 NHS ester at a molar ratio of 1:2 (protein:dye) by stirring in the dark.
5. Spin down the labeled nanoparticles by centrifugation at 8,000  $\times g$  for 10 min at room temperature, remove the supernatant and keep it to check labeling efficiency (step 3.7–3.8).
6. Wash the labeled nanoparticles by resuspending in 1 mL of 0.1 M sodium bicarbonate buffer by vortexing and pipetting up and down. Then discard the supernatant by centrifugation at 8,000  $\times g$  for 10 min at room temperature. Repeat step 3.6 for 5 times.
7. Resuspend the pellet in 5 mL of 0.1 M sodium bicarbonate buffer in a 15 mL centrifuge tube and make aliquots of the fluorescent nanoparticle in 1.5 mL centrifuge tubes (30 aliquots). Spin down at 8,000  $\times g$  for 10 min at room temperature, discard the supernatant and store at -80 °C protected from the light.
8. **Measure the labeling efficiency using a microvolume spectrophotometer.**
  1. Take 1  $\mu$ L of the original Atto 488 solution kept from step 3.3 and further dilute it in DMSO (e.g., 1:20 according to the volume ratio). This is the volume used to get the molar ratio of the protein nanoparticle and Atto 488 mix at step 3.4.
  2. Take 1  $\mu$ L of the saved supernatant from the labelling reaction in step 3.5. Measure the absorption (abs) at = 501 nm. The percentage of labelling is:
$$\left(1 - \frac{\text{abs supernatant}}{\text{abs Atto dilution}}\right) \times 100$$
9. Before the experiment, prepare the nanoparticle suspension at the desired concentration using 1x PBS solution.

## 4. Zebrafish Anesthetization and Oral Intubation

1. Fast the fish (>0.5 g) at least 48 h before the experiment to empty the intestine.
2. Move the fish (12 fish) to the experimental tanks (6 L) one night before the experiment to allow acclimatization<sup>12</sup>.
3. Vortex the nanoparticle solution well (e.g., 2,500 rpm and 30 s) and draw up the desired volume of nanoparticle suspension (e.g., 20–50  $\mu$ L) into the syringe attached to the protected needle.
4. Place the fish in the aerated 150 mg/L MS-222 solution (see section 2) until they sink to the bottom of the tank and do not respond to a tail fin pinch; this process takes less than 5 min.

5. Quickly transfer the anesthetized fish with a net to a wet plastic tray, orientate the animal horizontally to face the needle and immediately start the oral intubation.
6. Carefully support the fish with one hand and open the mouth with the other hand using the protected needle. Gently insert the needle down the esophagus to about 1 cm from the mouth opening.  
NOTE: The operator may feel a slight resistance when the end of the pipette tip has passed the gill. Take care not to angle the needle entry too much which may perforate the gill.
7. Slowly inject the nanoparticle suspension to the fish. Make sure the suspension does not flow outward through the gills or mouth.
8. Gently remove the needle and place the fish into the recovery tank (see section 2). Recovery usually takes within 1 min.
9. Check the fish carefully for any abnormality (e.g., bleeding at the gills is a sign of perforation).
10. Once the fish have recovered, return them to the experimental tanks.

## 5. Zebrafish Intestine Dissection

1. After a specified period of time post intubation (e.g., 5 h and/or 24 h), place the fish using a net into 300 mg/L MS-222 solution for euthanasia (see section 2). Make sure the operculum stops moving and there is no tail pinch reflex. Five minutes is normally sufficient.
2. Pick up the euthanized animal with a net and place it on a filter paper.  
NOTE: The filter paper is very useful for removing the adhesive tissue along the intestine.
3. Using sharp dissection scissors, make a semi-circular incision from the anus to the operculum and open incision using fine tweezers. Cut both ends of the intestine, take out all the internal organs, and place them on the filter paper.  
CAUTION: Work quickly to reduce cell metabolism and death.  
NOTE: Alternatively, remove the adhesive tissue in PBS and on ice.
4. Separate the intestine from internal organs making sure to keep its orientation (anterior to posterior intestinal segment) and stretch it out. Usually, the anterior segment of intestine is wider than posterior segment. Take care to obtain all of the intestine when dissecting.  
NOTE: The posterior end is quite fine and fragile in small fish and may break off, particularly in animals <0.7 g.
5. Roll the intestine on the filter paper with tweezers in order to detach the adhesive tissue from the intestine.
6. Proceed to prepare the intestine for various downstream analyses (sections 6, 7 and 8).

## 6. Preparing Intestinal Cells for Cytometry

1. Prepare in advance aliquots of 0.15% collagenase solution (see section 2.4).  
NOTE: Aliquots should be at room temperature before proceeding.
2. Optional: Continuing from step 5.5, slit open the intestine longitudinally and wash with 1x PBS.
3. Using tweezers, place the intestine in the 2 mL centrifuge tube filled with 0.15% collagenase solution.
4. Place the tubes on a vertical laboratory rotator for 1 h at room temperature in the dark.
5. Place the intestine on a 100  $\mu$ m cell strainer supported over a 50 mL centrifuge tube. Break up the intestine with a 5 mL syringe plunger, washing 3 times with 1x PBS, collecting the flow through sample in the 50 mL centrifuge tube.
6. Centrifuge the 50 mL centrifuge tube at 400  $\times$  g for 10 min, at 4 °C.
7. Carefully pipette off most of the supernatant while not losing the cells, some of which may be associated to the mucus.
8. Resuspend the intestinal cells at the bottom of the centrifuge tube with 500  $\mu$ L of 1x PBS and keep on ice until cytometry analysis.
9. Filter samples through a 30  $\mu$ m cell filter into 5 mL round bottom tubes for cytometry.
10. Set the parameters (e.g., the number of cells for analysis, the region of interest, voltage and compensation, selection of detectors) on a cytometer equipment (see materials).
11. Immediately analyze the cells on a cytometer, following the instructions of use<sup>13</sup>.

## 7. Preparing Intestine Cryosections for Confocal Microscopy

1. Fill the plastic mold (see **Table of Materials**) to half volume with optimum cutting temperature (O.C.T.) compound.
2. Continuing from step 5.5 and immediately after the dissection, carefully place the intestine in the plastic mold. Make sure the intestine is totally embedded in the O.C.T. compound. If necessary, add more O.C.T. compound to the plastic mold.  
NOTE: It is recommended to place the intestine with a "Z" shape in O.C.T. compound to easily follow its natural orientation.
3. Place the plastic mold on dry ice until it goes opaque (less than a minute).
4. Store the plastic mold at -80 °C for long-term use or process immediately using the following procedure.  
NOTE: The protocol can be paused here.
5. Slice the frozen intestine into 10  $\mu$ m sections or appropriate thickness using a cryostat at -20 °C.
6. Collect the intestine section with a fine brush onto a slide.
7. Immerse the slide in 4% PFA for 15 min at room temperature to fix the sample.  
CAUTION: PFA is toxic. Please read the material safety data sheet before working with it. Gloves and safety glasses should be worn, and always leave solutions inside a fume hood.
8. Wash the slide 3 times with 1x PBS, 10 min each.
9. Add a drop of mounting medium and place a coverslip over the specimen.  
NOTE: The protocol can be paused here.
10. Observe the sample under a confocal microscope at appropriate magnification.

## 8. Preparing the Intestine for Real Time qPCR (RT-qPCR)

1. Continuing from step 5.5, put the intestine in a cryogenic vial and rapidly freeze the intestine in liquid nitrogen and store at -80 °C until use.  
NOTE: The protocol can be paused here.

2. For homogenization, add 200  $\mu$ L of 2% (v/v) chilled 1-Thioglycerol/Homogenization solution (see **Table of Materials**) or alternative homogenization solution to the intestine sample.
3. Work quickly, homogenize the intestine sample on ice with a laboratory homogenizer at high speed (set at 25–30,000 rpm) until no visible tissue fragments remain. 3 times for 5 s is usually sufficient.
4. Isolate RNA using a commercial kit (see **Table of Materials**) according to the manufacturer's instructions<sup>14</sup> or a suitable alternative method. As needed, store the RNA at -80 °C for long-term use.  
NOTE: The protocol can be paused here.
5. Quantify the RNA concentration using a spectrophotometer<sup>15</sup> and evaluate the quality using a RNA analyzer<sup>16</sup>.
6. Prepare 1  $\mu$ g or an appropriate amount of cDNA using a cDNA synthesis kit according to the manufacturer's instructions.  
NOTE: The protocol can be paused here. For qPCR analysis, please follow the MIQE guidelines<sup>17</sup>.
7. Design appropriate primer pairs for the gene/s of interest.
8. Select a suitable reference gene and analyze the expression of each gene by a RT-qPCR detection system using a commercial kit (see materials).  
NOTE: For example, add 5  $\mu$ L of SYBR green supermix, 0.5  $\mu$ M primers, 2.5  $\mu$ L of diluted cDNA and 1.5  $\mu$ L of water in a final volume of 10  $\mu$ L for each well of qPCR plate.


## Representative Results

Zebrafish (average weight:  $1.03 \pm 0.16$  g) of mixed sex were successfully intubated with different recombinant protein nanoparticles (bacterial inclusion bodies) using our home-made oral intubation device (**Figure 1**). We have successfully performed the oral intubation and achieved a low average percentage mortality (6.8%) (**Table 1**). Zebrafish were either intubated with 30  $\mu$ L or 50  $\mu$ L of nanoparticle suspensions and the mortality was calculated within 24 h post intubation. The experiments were performed by two operators, R had less experience working with zebrafish oral intubation than J. The results showed that even a new operator could independently perform the oral intubation experiment and easily achieve a high survival rate by this protocol. From our experience, the optimal fish size is 1 g, but we have successfully intubated fish as small as 0.5 g.


To better understand whether the fluorescent nanoparticles IBs<sup>TNF $\alpha$</sup>  (a recombinant cytokine protein nanostructured as inclusion bodies) were delivered to zebrafish by our method and taken up by the zebrafish intestine or not, we performed cytometry analysis. The nanoparticles (100  $\mu$ g/fish, 50  $\mu$ L) and PBS (control) were orally administered (by intubation) to zebrafish and the intestine was dissected at 5 h and 24 h post intubation. Total intestinal cells were prepared by step 6 and analyzed by detecting the fluorescence emission signal. The representative histograms of fluorescence intensity and the dot plots of fluorescent cell percentage are shown in **Figure 2**. The density of fluorescent cells is clearly higher in nanoparticle intubated group compared to the control group in both 5 h and 24 h (**Figure 2A**). The percentages of fluorescent cells are significantly higher in both 5 h (46.3%) and 24 h (43.0%) of nanoparticle intubated groups (**Figure 2B**).

To further study which part of the intestinal layer is involved in the nanoparticle uptake, we performed confocal microscopy analysis. The nanoparticles (20  $\mu$ g/fish, 50  $\mu$ L) and PBS (control) were orally intubated to zebrafish and the intestine was dissected at 5 h post intubation. The intestine sections were prepared by a frozen tissue method according to step 7 (**Figure 3A**). The confocal images of fluorescent nanoparticles in the intestine are shown in **Figure 3B**. The fluorescent nanoparticles were found in zebrafish intestine. We observed the fluorescence in the epithelial cells, lamina propria, and muscle cells.


To verify whether we could extract high quality RNA by our protocol, we analyzed the RNA extracted from the intestine with a bioanalyzer (also referred to as RNA analyzer here). Zebrafish were intubated orally with PBS (50  $\mu$ L) or the nanoparticles (20  $\mu$ g/fish, 50  $\mu$ L). The intestines were dissected for RNA extraction at 24 h post intubation. We selected seven RNA samples to test with the analyzer. We found that all the tested RNA samples have high RNA integrity numbers ranging from 7.9 to 8.9 (**Figure 4**).



**Figure 1: Oral intubation device.** (A) Image of a 31 G Luer lock needle fixed on a 100  $\mu$ L syringe with the silicon tube and pipette tip end covering the tip of the needle. (B) An enlarged image of the needle part. The black arrow indicates where the pipette tip end exceeds the tip of the needle. [Please click here to view a larger version of this figure.](#)




**Figure 2: Flow cytometry analysis of fluorescent nanoparticles in zebrafish intestine via oral intubation.** Zebrafish were treated with PBS or fluorescent IBs<sup>TNFα</sup> (100 µg) for 5 h and 24 h, respectively. (A) Representative histograms of fluorescence intensity. (B) Dot plots graph of fluorescent cells percentage. Each green dot represents the percentage of fluorescent cells in one individual, n ≥ 4. Data represent mean ± standard error of the mean (SEM). Differences were analyzed using one-way ANOVA. Significant differences with respect to control (\*\*, p < 0.01) [Please click here to view a larger version of this figure.](#)



**Figure 3: Images of confocal microscopy analysis.** Fish were orally intubated with PBS or 20 µg/fish fluorescent nanoparticles. The intestine was dissected at 5 h post intubation. (A) Zebrafish intestine embedded in OCT compounds. The intestine was placed with the natural orientation of a "Z" shape (a: anterior end; p: posterior end). (B) Confocal microscopy images of zebrafish intestine. The white arrows show that the fluorescent nanoparticles are taken up in the intestinal mucosa. [Please click here to view a larger version of this figure.](#)





**Figure 4: RNA analyzer virtual gel image showing RNA extracted from 7 zebrafish intestines sampled post intubation.** Sample number 1 and 2 are PBS intubated groups and sample number 3 to 7 are nanoparticle intubated groups. RNA integrity numbers (RIN) are given at the bottom ranging from 7.9 to 8.9. [Please click here to view a larger version of this figure.](#)

Operator	Volume	# fish intubated	Mean weight (g ± SD)	# deaths	Mortality (%)
R	30 µL	22	0.88 ± 0.14	3	13.6
R	30 µL	17	0.93 ± 0.19	0	0
J	50 µL	19	1.23 ± 0.31	1	5.2
J	50 µL	30	1.08 ± 0.40	2	6.6
<b>Total</b>		<b>88</b>	<b>1.03 ± 0.16</b>	<b>6</b>	<b>6.8</b>

**SD: standard deviation of the mean**

**Table 1: Comparison of zebrafish mortality caused by two operators using the protocol.** The operator identification code, the intubated volume, the fish number, and the fish average weight in grams (g) are shown.

## Discussion

This protocol is an improvement of the previously described technique for oral intubation by Collymore *et al.*<sup>4</sup> Our protocol describes in detail the oral intubation method and includes the preparation of the intestine for downstream analyses. Our method improves fish manipulation speed allowing one person to perform the whole protocol rapidly, without much variation between operators. A main difference of our protocol with the previous one is that we evaluate the success of an oral intubation experiment not only by observation of the animal for well-being (e.g., no bleeding) and for no leakage of the fluid administered, but also by checking the uptake of a bioactive nanoparticle in the intestine using downstream analysis (cytometry, confocal microscopy and qPCR). We show that intubated fluorescently labeled nanoparticles were found in the intestine.

On the practical side, the intubation equipment is cheap, precise and reusable: the basic device is made of a reusable 100 µL glass (e.g., Hamilton) syringe coupled to a 31 G needle with a piece of silicon on the top. A cut sterile tip is placed over the silicon tube and can be changed for each individual administration. The syringe is reusable, and the thin needle allows intubating small fish with a high possibility of success. Moreover, the cut sterile tip could be placed directly over the needle without the silicon tube which allows the intubation equipment easier to be made in a laboratory. The administered bioactive compound is clearly visible, and a correct administration can be monitored easily. A critical step during the intubation procedure is the needle entry. It is important that the needle is not tilted or inserted too much to avoid perforating the gills. The mortality observed with this method is very low (approximately 7%) and depends on the size of the fish. Gill perforation is the most common cause of fish death and when it happens fish die within the first hour. Although fish of 0.5 g can be easily intubated the optimal fish size is around 1 g. Another difference with the method developed by Collymore *et al.* is that the fish are fasted for 48 h to be sure that gastrointestinal tract is empty. The whole intubation procedure can be done very quickly (30 animals/h) by one operator and importantly, the method is consistent between different operators<sup>4</sup>. The intubation method is easy to learn and does not require much practice to master.

The intestine dissection procedure must be properly performed to obtain good quality samples for cytometry, confocal and qPCR analyses. The critical step at this point is the dissection of the entire intestine; the posterior section is fragile and easy to lose. Once the intestine is dissected, it can be further processed for cytometry (2 h protocol), qPCR analysis (2 h protocol until total RNA isolation) or confocal microscopy (1 h to have samples ready for cryosection). It is very important to keep track of the orientation of the intestine, especially for confocal microscopy. Cytometry analysis requires rapid processing, and the procedure cannot be stopped before completion. Whereas, RNA isolation and samples prepared for cryosection can be properly stored and processed at any time. For cytometry, good quality samples without debris clumping the machine can be isolated with this method and fast analysis for the presence of fluorescent cells can be performed easily on individuals. RNA quality monitoring showed that high quality RNA can be isolated from intestine allowing analysis of individual fish by qPCR. Finally, the cryosection preparation for confocal microscopy provides important structural information about the protein nanoparticle uptake. Our methods thus provide a model to test the dynamics of oral probiotics or other compounds in the intestine.

The limitations of our study are the size of the fish since we did not test administration in fish smaller than 0.5 g and the use of a chemical anesthetic to sedate the animals. Some authors use cold water (0–4 °C) to anesthetize the zebrafish<sup>18</sup> but in the context of animal welfare and European legal constraints we decided that MS-222 was the method of choice.

The zebrafish offers many advantages over other model systems including a complete set of gene data and the available transgenic lines. For immunologists, the transgenic lines (e.g., *Tg mpx:GFP* and *Tg mpeg1:GFP*) set a scene for the *in vivo* observation of immune cells such as macrophages and neutrophils<sup>19,20</sup>. Combining our oral intubation and downstream analysis with transgenic lines, could be ideal for identifying the cell types involving the uptake, transportation and processing of oral vaccines, nanoparticles and pathogens in fish.

## Disclosures

The authors declare that no competing interests exist.

## Acknowledgements

This work was supported by grants from the Spanish Ministry of Science, European commission and AGAUR funds to NR (AGL2015-65129-R MINECO/FEDER and 2014SGR-345 AGAUR). RT holds a pre-doctoral scholarship from AGAUR (Spain), JJ was supported by a PhD fellowship from the China Scholarship Council (China) and NR is supported by the Ramón y Cajal program (RYC-2010-06210, 2010, MINECO). We thank Dr. Torrealba for expert advice in protein production, N. Barba from the "Servei de Microscopia" and Dr. M. Costa from the "Servei de Citometria" of the Universitat Autònoma de Barcelona for helpful technical assistance.

## References

1. Harrieff, M.J., Bermudez, L.E., Kent, M.L. Experimental exposure of zebrafish, *Danio rerio* (Hamilton), to *Mycobacterium marinum* and *Mycobacterium peregrinum* reveals the gastrointestinal tract as the primary route of infection: A potential model for environmental mycobacterial infection. *Journal of Fish Diseases*. **30** (10), 587-600 (2007).
2. Lovmo, S.D. *et al.* Translocation of nanoparticles and *Mycobacterium marinum* across the intestinal epithelium in zebrafish and the role of the mucosal immune system. *Developmental and Comparative Immunology*. **67**, 508-518 (2017).
3. Chen, T. *et al.* Small-Sized mPEG-PLGA Nanoparticles of Schisantherin A with Sustained Release for Enhanced Brain Uptake and Anti-Parkinsonian Activity. *ACS Applied Materials and Interfaces*. **9** (11), 9516-9527 (2017).
4. Collymore, C., Rasmussen, S., Tolwani, R.J. Gavage Adult Zebrafish. *Journal of Visualized Experiments*. (78), e50691-e50691 (2013).
5. Kim, S.H., Jang, Y.S. Antigen targeting to M cells for enhancing the efficacy of mucosal vaccines. *Experimental and Molecular Medicine*. **46** (3), e85 (2014).
6. Iwasaki, A., Medzhitov, R. Regulation of adaptive immunity by the innate immune system. *Science*. **327** (5963), 291-295 (2010).
7. Kunisawa, J., Kiyono, H. A marvel of mucosal T cells and secretory antibodies for the creation of first lines of defense. *Cellular and Molecular Life Sciences*. **62** (12), 1308-1321 (2005).
8. Rombout, J.H., Yang, G., Kiron, V. Adaptive immune responses at mucosal surfaces of teleost fish. *Fish Shellfish Immunology*. **40** (2), 634-643 (2014).
9. Salinas, I. The Mucosal Immune System of Teleost Fish. *Biology*. **4**, 525-39 (2015).
10. Munang'andu, H.M., Mutoloki, S., Evensen, Oø. An overview of challenges limiting the design of protective mucosal vaccines for finfish. *Frontiers in Immunology*. **6**, 542 (2015).
11. Lickwar, C.R. *et al.* Genomic dissection of conserved transcriptional regulation in intestinal epithelial cells. *PLoS Biology*. **15** (8), e2002054 (2017).
12. Brand, M., Granato, M., Nüsslein-Volhard, C. Keeping and raising zebrafish. *Zebrafish*. **261**, 7-37 (2002).
13. Rességuier, J. *et al.* Specific and efficient uptake of surfactant-free poly(lactic acid) nanovaccine vehicles by mucosal dendritic cells in adult zebrafish after bath immersion. *Frontiers in Immunology*. **8**, 190 (2017).
14. Kephart, D., Terry, G., Krueger, S., Hoffmann, K., Shenoi, H. High-Performance RNA Isolation Using the Maxwell 16 Total RNA Purification Kit. *Promega Notes*. (2006).
15. Thermo Fisher Scientific NanoDrop 1000 spectrophotometer V3.8 user's manual. *Thermo Fisher Scientific Incorporation*. (2010).
16. Lightfoot, S. Quantitation comparison of total RNA using the Agilent 2100 bioanalyzer, ribogreen analysis, and UV spectrometry. *Agilent Application Note*. (2002).
17. Huggett, J.F. *et al.* The digital MIQE guidelines: Minimum information for publication of quantitative digital PCR experiments. *Clinical Chemistry*. **59** (6), 892-902 (2013).
18. Matthews, M., Varga, Z.M. Anesthesia and euthanasia in zebrafish. *ILAR Journal*. **53** (2), 192-204 (2012).
19. Renshaw, S., Loynes, C. A transgenic zebrafish model of neutrophilic inflammation. *Blood*. **108** (13), 3976-3978 (2006).
20. Ellett, F., Pase, L., Hayman, J.W., Andrianopoulos, A., Lieschke, G.J. mpeg1 promoter transgenes direct macrophage-lineage expression in zebrafish. *Blood*. **117** (4) (2011).