

Materials List for:

The *Galleria mellonella* Waxworm Infection Model for Disseminated Candidiasis

Matthew J. Dunn¹, Andrew L. Woodruff¹, Matthew Z. Anderson^{1,2}

¹Department of Microbiology, The Ohio State University

²Department of Microbial Infection and Immunity, The Ohio State University

Correspondence to: Matthew Z. Anderson at anderson.3196@osu.edu

URL: <https://www.jove.com/video/58914>

DOI: [doi:10.3791/58914](https://doi.org/10.3791/58914)

Materials

Name	Company	Catalog Number	Comments
<i>Galleria mellonella</i>	Snackworms.com		Buy twice as many worms as expected to use
10 uL, Model 1701 N SYR Cemented needle, 26G, type 2 syringe	Hamilton	80000	
Petri dish, 100X15 mm, 500 pack	Fisher	FB0875712	
Microcentrifuge tube, 1.7 mL, 500 pack	VWR	87003-294	
Phosphate Buffered Saline (Biotechnology grade), 500 mL	VWR	97062-818	
Ethanol absolute, ≥99.5% pure, 500 mL	Millipore Sigma	EM-EX0276-1S	
autoclaved ddH ₂ O			