

Materials List for:

Isolating Malignant and Non-Malignant B Cells from *Ick:eGFP* Zebrafish

Jessica Burroughs-Garcia¹, Ameera Hasan^{*1}, Gilseung Park^{*1}, Chiara Borga¹, J. Kimble Frazer¹

¹Section of Pediatric Hematology-Oncology, Department of Pediatrics, University of Oklahoma Health Sciences Center

*These authors contributed equally

Correspondence to: J. Kimble Frazer at john-frazer@ouhsc.edu

URL: <https://www.jove.com/video/59191>

DOI: [doi:10.3791/59191](https://doi.org/10.3791/59191)

Materials

Name	Company	Catalog Number	Comments
35 µm mesh	Sefar Filter technology	7050-1220-000-13	
5 ml Polystyrene round-Bottom tube with cell-strainer cap	Falcon Corning Brand	352235	
50 ml conical tube	VWR international	525-0448	
AZ APO 100 Fluorescent microscope	Nikon		
Cytoflex	Beckman Coulter		
DS-Qi1MC camara	Nikon		
Ethyl 3-aminobenzoate methanesulfonate; MS-222	Sigma	E-10521	
FACSJazz	BD Biosciences		
Fetal bovine Serum	Thermo Fisher	10437028	
FlowJo v10.2	FlowJo, LLC		
<i>Ick:eGFP</i>			See Langeneu et al., 2004
NIS Elements software	Nikon	Version 4.13	
Penicilin -Streptomycin	Sigma	P4333	
Pestle micro-tube homogenizers	Electron Microscopy Sciences	64788-20	
Plastic Transfer pipettes			
<i>rag2:hMYC-ER</i>			See Gutierrez et al., 2011
RPMI Media 1640 1X	Life Technologies	11835-030	