

Materials List for:

Preparation of Neuronal Cultures from Midgastrula Stage Drosophila Embryos

Beatriz Sicaeros¹, Diane K. O'Dowd¹

¹Department of Development and Cell Biology, Department of Anatomy and Neurobiology, University of California, Irvine (UCI)

Correspondence to: Diane K. O'Dowd at dkodowd@uci.edu

URL: <http://www.jove.com/video/599>

DOI: [doi:10.3791/599](https://doi.org/10.3791/599)

Materials

Name	Type	Company	Catalog Number	Comments
Drosophila melanogaster	Animal			Fruit flies
Transferrin	Reagent	Sigma-Aldrich	T-1147	100x Stock: 10 mg/ml in water. Filter through 0.2 um syringe filter (cellulose acetate). Store in 220 ul aliquots (for 20 ml DMEM) at -20C.
Insulin		Sigma-Aldrich	I-6634	200x Stock: 10 mg/ml in 0.05N HCl. Filter through 0.2 um syringe filter (cellulose acetate). Store in 120 ul aliquots (for 20 ml DMEM) at -20C.
Putrescine		Sigma-Aldrich	P-5780	100x Stock: 10 mM in ddH ₂ O. Filter through 0.2 mm syringe filter (cellulose acetate). Store in 220 ml aliquots (for 20 ml DMEM) at -20°C
Selenium		Sigma-Aldrich	S-5261	100x Stock: 3 uM in ddH ₂ O. Put 0.0051 g Selenium in a 15 ml tube labeled A (3 mM stock). Add 10 ml of sterile water. Take 10 ul from Tube A to Tube B with 10 ml ddH ₂ O (3 uM stock). Filter Tube B through 0.2 um syringe filter (cellulose acetate). Store in 220 ul aliquots (for 20 ml DMEM) at -20C
Progesterone		Sigma-Aldrich	P-6149	100x Stock: 2 ug/ml1.Add 1ml of 100% EtOH to 0.001 g Progesterone bottle2.Add 49 ml ddH ₂ O3.Transfer 1 ml of this to second tube with 9 ml ddH ₂ O4.Filter through 0.2 um syringe filter (cellulose acetate)5.Store in 220 ul aliquots (for 20 ml DMEM) at -20C
DDM1	Medium			To 10 mls of DMEM add from stocks:100 ul Transferrin, 100 ul Putrescine, 100

				ul Selenium, 100 ul Progesterone, 50 ul Insulin
Petri dishes				
Cover-slips		Bellco Glass	1943-00012	Low lead glass, autoclaved.
Paper filter		Whatman, GE Healthcare	#1	
10cc syringe	Tool			
Agar				
Ethanol				
Acetic Acid				
Petri dishes				35 mm, plastic
Active baking yeast		Fleischman's		
10 cc syringe				
Buchner funnel				
Whatman #1 paper				
50% bleach				
Glass puller		Narishige International	PP-83	
Glass pipettes				
DMEM				Make DMEM: To 100 mls of Ham's F12/DME Media (DMEM) (Irvine Scientific #9052), add 0.476 g Hepes (20mM) (Sigma #H-3375). Mix. Add 1.25 ml L-Glutamine 200 mM (Irvine Scientific #9317). Mix. Filter in hood with 0.2 m acetate syringe filter (need 2 filters). The pH is 6.64 and Osm 288. Make aliquots of 10ml in 15ml centrifuge tubes. Store at 4°C for not more than 2 weeks.
The cultures made in this media must be maintained in a 4-5% CO2 incubator at about 22-24°C. We use a standard mammalian tissue culture incubator placed in a cold room and heated to 23°C. Always used autoclaved filtered water and make in containers that are reserved for media and supplements only. Never put a pH meter or stir bar used for other purposes in media.				