

Video Article

# Monitoring Actin Disassembly with Time-lapse Microscopy

Hao Yuan Kueh<sup>1</sup>

<sup>1</sup>Dept. of Systems Biology, Harvard Medical School

Correspondence to: Hao Yuan Kueh at [kueh@fas.harvard.edu](mailto:kueh@fas.harvard.edu)

URL: <http://www.jove.com/video/66>

DOI: [doi:10.3791/66](https://doi.org/10.3791/66)

Keywords: Cellular Biology, Issue 1, cytoskeleton, actin, timelapse, filament, chamber

Date Published: 11/8/2006

Citation: Kueh, H.Y. Monitoring Actin Disassembly with Time-lapse Microscopy. *J. Vis. Exp.* (1), e66, doi:10.3791/66 (2006).

## Abstract

## Video Link

The video component of this article can be found at <http://www.jove.com/video/66/>

## Protocol

1. Construct a flow chamber as shown in the video. Make sure that the parafilm seal is tight and use washed coverslips.
2. Incubate actin-binding agent in the chamber for 5-10'.
3. Block non-specific binding sites on the glass coverslip with a blocking protein.
4. Polymerize actin inside the chamber by flowing in G-actin in polymerizing buffer.
5. Washout unpolymerized actin by flowing in excess buffer.