

Oligo name	Length in nt	Sequence	Purpose	Manufacturer (Purification)
full Kwok + MiSeq	61	5'-PHO-tgaagagcctagtcgctgttcannnnnnctgccatagagagatcggaagagcacacgtct-SpC3-3'	Adaptor	IDT DNA technologies (HPLC)
2xBiotin SS bait	40	5'-biotin-cagtgtggaaaatctctagcagtggcgcccgaacagggac-biotin-3'	Hybrid Capture	MWG Eurofins (HPLC)
Biotin 1-16 ss	22	5'-cagtgtggaaaatctctagcag-BiTEG-3'	Hybrid Capture	MWG Eurofins (HPLC)
Biotin tRNA + CTG	16	5'-cagtggcgcccgaaca-BITEG-3'	Hybrid Capture	MWG Eurofins (HPLC)
MP1.0 + 22HIV	82	5'-aatgatacggcgaccaccgagatctacactctttccctacacgacgctcttccgatctcactgctagagattttccacactg-3'	PCR amplification	MWG Eurofins (HPLC)