

Current Stage Position (mm) X: 8.086 Y: 22.026 Jog Stage

Slide 131 [TMA569-MD]



Current stage position: Target stage position:

Target Point:

Target Selection:

Stage to Target

Jog Stage

Jog Direction

<<<

<<

<

>

>>

>>>

<<<

<<

<

>

>>

>>>

>>>

>>

>

<

<<

<<<

>>>

>>

>

<

<<

<<<

Jog Increment Distance (μm)

10

Close

Mode: SED

Imaging Mode: QC - 300μm

Reapply

Reset

Sample Current (nA): -0.3519

Gain: 57.769 set

Coarse:

Fine:

SED Image

Frame Size: 512x512

Dwell (μs): 10

FOV Size (μm): 1800

50 1850

Focus

Current Focus (kV): 20.053 set

Coarse:

Medium:

Fine:

Save

Stigmation

Set Levels

Fit to window (100%)

Fit to window

Crosshair

Expand

FOV Center Point:

X: 8.086

Y: 22.026

Add FOV

Download SED

Estimated Run Time (hh:mm:ss):

Sample Name

FOV Center Point

Scans

Frame

Timing Choice

Imaging Mode

FOV Size (μm)

Section / Target

Estimate

Contact support