

Video Article

Isolation of Primary Myofibroblasts from Mouse and Human Colon Tissue

Hassan Khalil¹, Wenxian Nie¹, Robert A Edwards², James Yoo¹

¹Department of Surgery, UCLA

²Department of Pathology, UC Irvine

Correspondence to: James Yoo at JaYoo@mednet.ucla.edu

URL: <https://www.jove.com/video/50611>

DOI: [doi:10.3791/50611](https://doi.org/10.3791/50611)

Keywords: Cellular Biology, Issue 80, Myofibroblasts, Mesenchymal Stromal Cells, Gastrointestinal Tract, stroma, colon, primary cells

Date Published: 10/12/2013

Citation: Khalil, H., Nie, W., Edwards, R.A., Yoo, J. Isolation of Primary Myofibroblasts from Mouse and Human Colon Tissue. *J. Vis. Exp.* (80), e50611, doi:10.3791/50611 (2013).

Abstract

The myofibroblast is a stromal cell of the gastrointestinal (GI) tract that has been gaining considerable attention for its critical role in many GI functions. While several myofibroblast cell lines are commercially available to study these cells *in vitro*, research results from a cell line exposed to experimental cell culture conditions have inherent limitations due to the overly reductionist nature of the work. Use of primary myofibroblasts offers a great advantage in terms of confirming experimental findings identified in a cell line. Isolation of primary myofibroblasts from an animal model allows for the study of myofibroblasts under conditions that more closely mimic the disease state being studied. Isolation of primary myofibroblasts from human colon tissue provides arguably the most relevant experimental data, since the cells come directly from patients with the underlying disease. We describe a well-established technique that can be utilized to isolate primary myofibroblasts from both mouse and human colon tissue. These isolated cells have been characterized to be alpha-smooth muscle actin and vimentin-positive, and desmin-negative, consistent with subepithelial intestinal myofibroblasts. Primary myofibroblast cells can be grown in cell culture and used for experimental purposes over a limited number of passages.

Video Link

The video component of this article can be found at <https://www.jove.com/video/50611/>

Introduction

The regulation of gastrointestinal (GI) tract function involves complex, dynamic interactions between the mucosal layer and its underlying mesenchyme. These interactions normally serve to maintain homeostasis but may also lead to the development of disease states under pathologic conditions¹. Myofibroblasts are a subpopulation of stromal cells located subjacent to the epithelial layer that communicate in a paracrine fashion with surrounding cell subpopulations to regulate a number of critical GI tract processes that include mucosal regeneration, repair, and fibrosis².

The 18Co cell line is an example of a human colonic myofibroblast cell line that has been widely used to study myofibroblast function because it shares many of the characteristics of primary *in situ* myofibroblasts. These include the protein expression of α -smooth muscle actin (α -SMA) and vimentin, as well as the expression of a number of cell surface receptors such as the epidermal growth factor receptor, or receptors for lysophosphatidic acid^{3,4}. Because of these shared ultrastructural characteristics, as well as similarities in biologic function⁵, the 18Co cell line has been used extensively to study myofibroblast function in the context of many disease states such as inflammatory bowel disease or colorectal cancer^{3,6}. However, there are inherent limitations and concerns when using a cell line. These include genotypic instability over time that differentiates cell lines from primary cells, altering the phenotype and biologic function of the cell, including growth rates and interactions with other cell populations. Cell lines also lack the normal components of the GI microenvironment (epithelial, stromal, and vascular components), which is a major limitation to its use. Therefore, the conclusions drawn from experimental cell line research requires further validation to reduce the risk of misinterpretation.

Primary myofibroblasts can be obtained from human and mouse colon tissue to confirm experimental findings identified in a cell line⁷. The method was originally described by Mahida, *et al.*⁸, and our protocol is one of many well-described techniques that can be used to isolate primary myofibroblasts from mouse and human colon⁹. For any mouse model of colonic disease, this technique can be used to isolate primary myofibroblasts to study how they interact with neighboring cells or contribute to both normal or pathologic GI processes^{10,11}. This technique can be used to study how myofibroblasts contribute to mucosal healing, fibrosis, and the development of colorectal adenomas and carcinomas. Primary myofibroblasts can also be obtained from human colon tissue that has been resected at the time of operation for benign (inflammatory bowel disease, strictures, diverticulitis) or malignant conditions. Isolated primary cells from human and mouse colon tissue can be grown in cell culture and utilized over a limited number of passages.

Protocol

1. Obtain Colon Tissue

Mouse

A protocol to perform animal experiments, detailing the rationale and objectives of the research, was submitted and approved by our institutional Office of Animal Research Oversight.

1. Euthanize the mouse with inhaled isoflurane followed by cervical dislocation.
2. Make a ventral midline incision, retract the small bowel away and identify the colon. Retract the bladder and uterus (if present) laterally and identify the pubic bone and cut it with fine scissors to expose the rectum.
3. Dissect the dorsal mesentery of the rectum off of the colon and mobilize the colon from the anus to the cecum, and transect it to disconnect the colon from the small bowel.
4. Place the specimen in ice-cold phosphate buffered saline (PBS).
5. Fill a 10 ml syringe with ice-cold PBS and flush the colon lumen several times to clear its contents, then cut the colon open length-wise using scissors.
6. Place the colon in a 50 ml conical tube containing 20 ml of ice-cold PBS. The colon does not have to be placed in any particular orientation. Wash the colon by replacing the PBS until the fluid is clear with no particulate matter.

Human

A protocol to obtain human tissue from surgical patients, detailing the importance of obtaining this tissue as well as an assessment of potential risks and benefits, was submitted and approved by our institutional Office of the Human Research Protection Program.

1. Research collaboration with a colorectal surgeon is necessary to be able to obtain human tissue for your research. A protocol to obtain human colon tissue should be submitted and approved by your institutional review board.
2. Written informed consent is obtained from the patient prior to the surgical procedure.
3. The surgeon will perform a colon resection in the standard fashion. Once the colon has been removed from the patient, the specimen is immediately taken to pathology. A 0.5 in full-thickness section of colon wall not needed for the pathologic evaluation is provided and is immediately placed on ice-cold PBS. This specimen should not include the colon mesentery, and epiploic appendages should be removed as well, since the presence of fat will affect the digestion process.

2. Denude Epithelial Cells

1. Incubate the entire mouse colon in 25 ml of 5 mM EDTA in HBSS (room temperature) in a 50 ml conical tube. Place the conical tube in a 37 °C shaking air bath (250 rpm) for 15 min. After 15 min, carefully pour off the fluid and refill the tube with 25 ml of 5 mM EDTA, repeat for a total of 5 washes. For human colon, after obtaining a ½ inch strip of colon wall from the pathologist cut the colon with scissors into four equal-sized pieces. Place two of the pieces into one 50 ml conical tube, and place the remaining two pieces of colon tissue into a separate 50 ml conical tube. Then incubate with 25 ml of 5 mM EDTA in HBSS and perform washes, as described above.
2. Rinse the colon twice in 20 ml of ice-cold PBS.
3. Place the mouse colon tissue in a new 50 ml conical tube containing 20 ml of RPMI-5, 10 U of dispase, and 2,000 U of collagenase D (room temperature). For the human colon tissue, use twice the amount of dispase and collagenase (20 U dispase, 4,000 U collagenase D).
4. Place the tube in a 37 °C shaking air bath oriented vertically (If placed horizontally, there is too much aeration and cell damage). Shake at 250 rpm for up to 60 min. After exposure to the dispase and collagenase D, the colon tissue will start to digest and the colon tissue will begin to look stringy. This usually takes around 30 min. When the colon has this appearance, remove the colon from the enzymes.
5. Shake each tube up and down 2-4 times to break up the tissue. Pellet the tissue at 200 x g in a tabletop centrifuge (4 °C) for 5 min. Carefully pour off the supernatant and discard.

3. Myofibroblast Isolation and Culture

1. Resuspend each pellet by adding 10 ml ACK lysis buffer (47 °C). Centrifuge again at 200 x g for 5 min (4 °C), pour off and discard the supernatant.
2. Resuspend each pellet by adding 10 ml RPMI-5 with a pipette.
3. Pass half of the cell suspension (5 ml) through a 70 µm mesh strainer using a 10 ml pipette into a 100 mm TC-treated dish. Then add an additional 15 ml of RPMI-5. Repeat with the remaining 5 ml of cell suspension into a different 100 mm TC-treated dish. One mouse colon will therefore yield two 100 mm TC-treated dishes. One ½ inch strip of human colon tissue will yield a total of four 100-mm TC-treated dishes.
4. Place the tissue culture dishes in a 10% CO₂ incubator at 37 °C. The culture conditions are the same for mouse and human primary cells.
5. After 3 hr, gently wash off non-adherent cells with two changes of HBSS. The debris floating should be gently washed off. Adherent cells are composed of epithelial cells, macrophages and myofibroblasts. One week post-seeding, only myofibroblasts divide, macrophages and epithelial cells senesce after the first passage. Add fresh RPMI-5 and grow the cells in cell culture.
6. Cells are passaged by trypsinization for 5 min and are split 2:1.

Abbreviations

HBSS: Hank's balanced salt solution; RPMI-5: Roswell Park Memorial Institute media; ACK: Ammonium-Chloride-Potassium; FCS: fetal calf serum; TC-treated: tissue culture treated


Solutions

ACK lysis buffer - (4.15 g NH_4Cl , 0.5 g KHCO_3 , 18.6 mg Na_2EDTA , 400 ml H_2O , adjust pH to 7.2-7.4 with 1 N HCl). Filter sterilize through a 0.2 μm filter and store at 4 °C. **RPMI-5** - (to make 500 ml: 454.5 ml RPMI, 25 ml FCS, 5 ml 200 mM L-glutamine, 5 ml 1 M HEPES, pH 7.4, 5 ml 100 mM sodium pyruvate, 5 ml 100x Pen-Strep, 500 ml 50 mM B-ME in PBS).


Representative Results

Once isolated, primary human myofibroblasts can be grown in cell culture and used over a limited number of passages (up to passage 4). These cells are characterized as being α -smooth muscle actin and vimentin positive, and desmin-negative (**Figure 1A**), consistent with intestinal subepithelial myofibroblasts^{5,7}. They also have a characteristic stellate morphology (**Figure 1B**). Primary myofibroblasts can then be used to confirm experimental findings identified in a cell line. This approach was utilized to demonstrate that the pro-inflammatory cytokine TNF- α (10 ng/ml) induced the upregulation of EGF receptor expression and signaling in these cells⁷. Upregulated EGF receptor expression and signaling was initially found utilizing the previously mentioned human colonic myofibroblast cell line 18Co (**Figure 2**). In **Figure 2A**, we demonstrate that exposure of 18Co cells to TNF- α led to a time-dependent increase in EGF receptor expression. This correlated with enhanced EGF-induced COX-2 expression and p42/44 MAPK phosphorylation in these cells (**Figure 2B**). To validate these experimental findings, experiments were repeated using primary human myofibroblasts isolated from surgically resected colon of patients with colorectal cancer. As shown in **Figure 2C** and **2D**, primary human myofibroblasts closely mimic the experimental findings initially identified in the 18Co cell line⁷.

1A




← α -SMA (42-kDa)



← Vimentin (57-kDa)

1B



Ab:

Figure 1. Primary myofibroblasts can be isolated from human colon tissue⁷. Isolated primary cells demonstrate a myofibroblast-like phenotype that are consistently α -SMA and vimentin positive (**Figure 1A**) and demonstrate a stellate morphology (**Figure 1B**).




Figure 2. Experimental findings in a cell line are substantiated using primary human myofibroblast cells⁷. The human colonic myofibroblast cell line 18Co was used to demonstrate that TNF- α induces the upregulation of EGF receptor expression and signaling in these cells (Figure 2A). This upregulation of EGF receptor expression was associated with enhanced EGF receptor signaling, shown in Figure 2B, where subsequent exposure to EGF led to enhanced p42/44 MAPK phosphorylation and COX-2 expression. These effects were completely inhibited with the EGF receptor inhibitor AG1478. Primary human myofibroblasts isolated from surgically resected colon of patients with colorectal cancer confirm these experimental findings (Figures 2C and 2D). Tumor necrosis factor- α =TNF- α , EGFR=EGF receptor, AG=AG1478, COX-2=cyclo-oxygenase-2. [Click here to view larger figure.](#)




Figure 3. Primary myofibroblasts can be isolated from mouse colon tissue. A 10X view of primary mouse myofibroblast cells. [Click here to view larger figure.](#)

Discussion

The general technique is similar when using either mouse or human colon. This technique can also be used to isolate myofibroblasts from the small intestine. Specific details for the isolation of myofibroblasts from 1. mouse and 2. human colon are discussed below.

1. Mouse Colon

Irrigation of the mouse colon is facilitated by attaching a 4 in, 16 G Popper blunt-end needle to one end of the colon. This also helps when cutting the colon lengthwise. You can accordion the colon tissue onto the needle, then place one tip of the scissor with the sharp end up and slide the colon tissue over the scissor edge while cutting to open it lengthwise.

In step 2.1, it is important to understand that nuclear material from dead cells can lead to cell clumping. The addition of DNase (50 mg/ml) may be used to address this problem.

The most critical step is step 2.4. Keep a careful watch of the colon tissue. When the colon begins to look stringy, remove the colon from the enzymes. You do not need to incubate the colon for the entire 60 min. If the colon tissue is exposed to the dispase and collagenase D for too long, the cells will die. This is a common mistake that will result in a low cell yield.

Contamination of the cells is another common problem. Careful technique, appropriate washing and filtering will minimize the risk of contamination, along with standard sterile techniques used for cell culture work.

2. Human Colon

The length of time for incubation with the enzymes depends on the size of the specimen that is being digested. We typically receive a ½ inch full thickness strip of human colon wall. There is a modification in step 2.3 when isolating myofibroblasts from human colon tissue. For step 2.3, for human colon you must double the amount of dispase (20 U) and collagenase D (4000 U) used. An even larger segment of colon will require a longer incubation time. Also, keep in mind that the specific activity of the enzymes may vary with lots, so you may have to adjust the incubation time accordingly.

For both primary mouse and human myofibroblasts, Once grown in culture, the isolated cells should be evaluated to assess for purity of the cells and to confirm that the cells are myofibroblast-like. This can be done with a variety of technique including immunohistochemical staining (antibodies to myofibroblast markers include α -smooth muscle actin and vimentin; monoclonal antibody to CD68 is specific to macrophages, CD31 is specific to endothelial cells, and CD45 is specific to immune cells; antibodies to e-cadherin and various cytokeratins are specific for epithelial cells. mRNA expression can also be assessed by RT-PCR¹². Flow cytometry can also be used for cell analysis¹³.

Primary myofibroblasts can be studied, frozen in liquid nitrogen for subsequent use at a later date, or used for *in vitro* and co-culture experimentation. Primary cells can also be implanted into the subcutaneous tissue of mice to study cell-cell interactions¹⁴. Future applications may involve re-implantation of primary myofibroblasts into the colon wall of a syngeneic mouse utilizing mouse colonoscopy following genetic manipulation.

A number of alternative methods to isolate primary myofibroblasts from human and mouse colon have also been described. Two are listed below:

1. Mahida, *et al. Am. J. Physiol.* **273** G1341-1348, (1997).
2. De Wever, *et al. Int. J. Oncol.* **39**, 393-400, (2011).

Disclosures

Authors have nothing to disclose.

Acknowledgements

This work was supported by National Institutes of Health Grant 5KO8DK085136-02 to JY.

References

1. Armaka, M. *et al.* Mesenchymal cell targeting by TNF as a common pathogenic principle in chronic inflammatory joint and intestinal diseases. *J. Exp. Med.* **205**, 331-337 (2008).
2. Powell, D. W. *et al.* Myofibroblasts. I. Paracrine cells important in health and disease. *Am. J. Physiol.* **277**, C1-9 (1999).
3. Yoo, J., Rodriguez Perez, C. E., Nie, W., Sinnott-Smith, J. & Rozengurt, E. Protein kinase D1 mediates synergistic MMP-3 expression induced by TNF- α and bradykinin in human colonic myofibroblasts. *Biochem. Biophys. Res. Commun.* **413**, 30-35 (2011).
4. Yoo, J., Chung, C., Slice, L., Sinnott-Smith, J. & Rozengurt, E. Protein kinase D mediates synergistic expression of COX-2 induced by TNF- α and bradykinin in human colonic myofibroblasts. *Am. J. Physiol. Cell Physiol.* **297**, C1576-1587, doi:10.1152/ajpcell.00184.2009 (2009).
5. Valentich, J. D., Popov, V., Saada, J. I. & Powell, D. W. Phenotypic characterization of an intestinal subepithelial myofibroblast cell line. *Am. J. Physiol.* **272**, C1513-1524 (1997).

6. Rodriguez Perez, C. E., Nie, W., Sinnott-Smith, J., Rozengurt, E. & Yoo, J. TNF-alpha potentiates lysophosphatidic acid-induced COX-2 expression via PKD in human colonic myofibroblasts. *Am. J. Physiol. Gastrointest. Liver Physiol.* **300**, G637-646 (2011).
7. Yoo, J. *et al.* TNF-alpha induces upregulation of EGFR expression and signaling in human colonic myofibroblasts. *Am. J. Physiol. Gastrointest. Liver Physiol.* **302**, G805-814 (2012).
8. Mahida, Y. R. *et al.* Adult human colonic subepithelial myofibroblasts express extracellular matrix proteins and cyclooxygenase-1 and -2. *Am. J. Physiol.* **273**, G1341-1348 (1997).
9. De Wever, O. *et al.* Priming and potentiation of DNA damage response by fibronectin in human colon cancer cells and tumor-derived myofibroblasts. *Int. J. Oncol.* **39**, 393-400, doi:10.3892/ijo.2011.1034 (2011).
10. Edwards, R. A. & Smock, A. Z. Defective arachidonate release and PGE2 production in Gi alpha2-deficient intestinal and colonic subepithelial myofibroblasts. *Inflamm. Bowel Dis.* **12**, 153-165 (2006).
11. Hoang, B., Trinh, A., Birnbaumer, L. & Edwards, R. A. Decreased MAPK- and PGE2-dependent IL-11 production in Gialpha2-/- colonic myofibroblasts. *Am. J. Physiol. Gastrointest. Liver Physiol.* **292**, G1511-1519 (2007).
12. Li, Z. B., Kollias, H. D. & Wagner, K. R. Myostatin directly regulates skeletal muscle fibrosis. *J. Biol. Chem.* **283**, 19371-19378, doi:10.1074/jbc.M802585200 (2008).
13. Saada, J. I. *et al.* Subepithelial myofibroblasts are novel nonprofessional APCs in the human colonic mucosa. *J. Immunol.* **177**, 5968-5979 (2006).
14. Lahar, N. *et al.* Intestinal subepithelial myofibroblasts support *in vitro* and *in vivo* growth of human small intestinal epithelium. *PLoS One.* **6**, e26898 (2011).