Materials List for:

MRI-guided Focused Ultrasound Thalamotomy for Patients with Medically-refractory Essential Tremor

Ying Meng¹, Yuexi Huang², Benjamin Solomon², Kullervo Hynynen², Nadia Scantlebury¹, Michael L. Schwartz¹, Nir Lipsman¹

Correspondence to: Nir Lipsman at nir.lipsman@sunnybrook.ca

URL: https://www.jove.com/video/56365

DOI: doi:10.3791/56365

Materials

Name	Company	Catalog Number	Comments
InSightec ExAblate Neuro system	InSightec		The ExAblate Neuro consists of a helmet transducer with 1,024 independent elements, each with central frequency 650 kHz and is integrated with a 3 Tesla MR scanner.
3 Tesla MRI scanner			
MRI compatible Cosman-Roberts- Wells (CRW) stereotactic frame			
20 mL of 1% lidocaine with 1:100,000 epinephrine			
hair shaver/clipper			
warming blanket			
compression stockings			

¹Division of Neurosurgery, Sunnybrook Health Sciences Centre

²Sunnybrook Research Institute, Sunnybrook Health Sciences Centre